


## The Chamberlain

A Vehicle Seat Belt positioning device designed for Children and Adults 81-225 lbs. with special health care needs

Keeps the occupant in position while using the vehicle seat belt

## Design Features

- Head support
- Small Medium and large vests for torso support
- Special vest sizes available
- Vest attached to harness for easy installation
- Guide to keep vehicle belt in position on shoulder and off the neck



**The Chamberlain**

## The Chamberlain



**Harness- vest and headrest mounts to vehicle seat**

# The Chamberlain



**Harness and Vest remains in vehicle when not in use**

# Design Features



Leg-Hip straps



Prevents slouching under the vehicle seat belt

## The Chamberlain

# Design Features

- E-Z –UP Cap available on the Chamberlain
- Cap attaches to Velcro located on Chamberlain headrest to keep head from falling forward




- Three cap sizes available
- Caps adjustable to fit child's head

## The Chamberlain

## Design Features

- No base required
- Occupant sits on vehicle seat
- Easy transfer into vehicle
- Children and Adults  
81 to 225 lbs.



**The Chamberlain**