



2022 Chamberlain Retail Price List

The Chamberlain-Pediatric and Adult Positioning Vehicle System and Accessories

The Chamberlain is available in three fabric color options. To select a color, simply replace the asterisk in the part numbers in the table with the color option abbreviation from the following list.

- Sapphire Blue – S
- Total Camo – TC
- Beige Naugahyde - N

See Important Statement Below About Installing the Chamberlain***

Part Number	Description	Retail Price
4000CHS-*	The Chamberlain Standard* – an adaptive pediatric and adult positioning vehicle system (Weight range: 81 to 225 lbs. A Positioning Vest is included in price. Please choose size below.	\$795.00
2000CVS-*	Small Positioning Vest – for chest circumference between 25 and 32 inches and belly circumference less than 32 inches. **Price charged only when ordering vest separately.	\$375.00**
2000CVM-*	Medium Positioning Vest – for chest circumference between 30 and 37 inches and belly circumference less than 37 inches. **Price charged only when ordering vest separately.	\$375.00**
2000CVL-*	Large Positioning Vest – for chest circumference between 35 and 44 inches and belly circumference less than 44 inches. **Price charged only when ordering vest separately.	\$375.00**
2000VB	Add Side Release Buckles to one of the vests above. (Side release buckles can help stop an occupant from escaping the vest).	\$99.00
	Optional Accessories	
2000EUHR-*	EZ-Up Head Rest - Supports head to limit side to side motion. To limit forward motion, choose a cap size from below:	195.00
See note to the right regarding the three cap material options.	Note: There are three options for cap material. Sapphire Blue, Total Camo, or Khaki. The Khaki option is the coolest option for hot temperatures because it is single layer.	
1000EUCS-K	EZ-Up Cap – Size Small (Khaki) - (Head Circumference: 17.5 – 19.5 inches)	\$99.00
1000EUCS-*	EZ-Up Cap – Size Small - (Head Circumference: 17.5 – 19.5 inches)	\$99.00
1000EUCM-K	EZ-Up Cap – Size Medium (Khaki) - (Head Circumference: 17.5 – 19.5 inches)	\$99.00

1000EUCM-*	EZ-Up Cap – Size Medium - (Head Circumference: 19.5 – 21.5 inches)	\$99.00
1000EUCL-K	EZ-Up Cap – Size Large (Khaki) - (Head Circumference: 17.5 – 19.5 inches)	\$99.00
1000EUCL-*	EZ-Up Cap – Size Large - (Head Circumference: 17.5 – 19.5 inches)	\$99.00
2000AB	Attachment Bracket Kit (for Vehicles without LATCH) – See Important Note Below ***	\$95.00

* The “Chamberlain Standard” includes the following features:

- Torso Support System
- Anti-Slouching Leg Strap Assembly
- Positioning Vest (Small, Medium, or Large)
- Tether Extension Assembly

**One Vest is included in price of “Chamberlain Standard with Vest.” Retail Price of additional vests is \$375.00 each.

*** The Attachment Bracket Kit is designed for use in vehicles that do not have LATCH. The Chamberlain must be installed with LATCH or with the Attachment Bracket kit. All Vehicles 2003 and newer are equipped with LATCH in at least two rear seating positions, so there should be no installation issues for these vehicles and the Attachment Bracket will not be necessary. Many 2000 - 2002 models also have LATCH (See your vehicle owner’s manual to verify if LATCH is present).

For vehicles without LATCH, the attachment bracket can only be used on vehicle seats that you can slide a ruler, or similar object, from the back side of the vehicle seat through the crack in the vehicle seat (where the seat bottom and seat back come together) through to the front side of the seat. This can usually be done in mini-vans and SUVs. Generally, the Attachment Bracket cannot be installed in the back seat of sedans. Installation of the Attachment Bracket Kit does require that a hole be drilled through the vehicle’s floor board.

Shipping charges – The cost of shipping “Best Way” will be added to Purchase Orders.
--

Thanks for Your Business!