



2022 Churchill Retail Price List

The Churchill-Pediatric Positioning Vehicle Restraint System and Accessories

See Important Statement Below About Installing the Churchill****

The Churchill is available in three fabric color options. To select a color, simply replace the asterisk in the part numbers in the table with the color option abbreviation from the following list.

- Sapphire Blue – S
- Total Camo – TC
- Beige Naugahyde - N

Part Number	Description	Retail Price
2000CS15-*	The Churchill Standard* with Vest and 15 Inch Seat Depth – a pediatric positioning vehicle restraint system (Weight range: 44 to 175 lbs.; Height Range: 44 to 72 inches). One Positioning Vest is included in price. Please choose vest size below.	\$1,095.00
2000CS18-*	The Churchill Standard* with Vest and 18 Inch Seat Depth – a pediatric positioning vehicle restraint system (Weight range: 44 to 175 lbs.; Height Range: 44 to 72 inches). One Positioning Vest is included in price. Please choose vest size below.	\$1,095.00
2000CH15-*	The Churchill with Positioning Harness*** and 15 Inch Seat Depth - <u>CAN BE USED WITH OUR ESCAPE PROOF PRODUCTS</u> – a pediatric positioning vehicle restraint system (Weight range: 44 to 175 lbs.; Height Range: 44 to 72 inches). One Positioning Harness is included in price. NOTE: Add the escape proof Buckle Guard and Chest Clip Guard listed below for the proficient escape artist.	\$995.00
2000CH18-*	The Churchill with Positioning Harness*** and 18 Inch Seat Depth - <u>CAN BE USED WITH OUR ESCAPE PROOF PRODUCTS</u> – a pediatric positioning vehicle restraint system (Weight range: 44 to 175 lbs.; Height Range: 44 to 72 inches). One Positioning Harness is included in price. NOTE: Add the escape proof Buckle Guard and Chest Clip Guard listed below for the proficient escape artist.	\$995.00
2000CVS-*	Small Positioning Vest - for chest circumference between 25 and 32 inches and belly circumference less than 32 inches. **Price charged only when ordering vest separately.	\$375.00**
2000CVM-*	Medium Positioning Vest - for chest circumference between 30 and 37 inches and belly circumference less than 37 inches. **Price charged only when ordering vest separately.	\$375.00**
2000CVL-*	Large Positioning Vest - for chest circumference between 35 and 44 inches and belly circumference less than 44 inches. **Price charged only when ordering vest separately.	\$375.00**

2000VB	Add Side Release Buckles to one of the vests above. (Side release buckles can help stop an occupant from escaping the vest. For the proficient escape artist, rather than the vest with buckles, we suggest selecting the harness version with both escape products below).	\$99.00
	Optional Accessories	
2000BG	Buckle Guard – ONLY FOR USE WITH 2000CH15-* or 2000CH18-*	\$150.00
2000CCG-*	Chest Clip Guard – ONLY FOR USE WITH 2000CH15-* or 2000CH18-*	\$150.00
2000EUHR-*	EZ-Up Head Rest – Supports head to limit side to side motion. To limit forward motion, choose a cap size from below: The cap is designed to Velcro to the head rest.	\$195.00
See note to the right regarding the three cap material options.	Note: There are three options for cap material. Sapphire Blue, Total Camo, or Khaki. The Khaki option is the coolest option for hot temperatures because it is single layer.	
1000EUCS-K	EZ-Up Cap – Size Small (Khaki) - (Head Circumference: 17.5 – 19.5 inches)	\$99.00
1000EUCS-*	EZ-Up Cap – Size Small - (Head Circumference: 17.5 – 19.5 inches)	\$99.00
1000EUCM-K	EZ-Up Cap – Size Medium (Khaki) - (Head Circumference: 19.5 – 21.5 inches)	\$99.00
1000EUCM-*	EZ-Up Cap – Size Medium (Sapphire Blue) - (Head Circumference: 19.5 – 21.5 inches)	\$99.00
1000EUCL-K	EZ-Up Cap – Size Large (Khaki) - (Head Circumference: 21 – 23.25 inches)	\$99.00
1000EUCL-*	EZ-Up Cap – Size Large (Sapphire Blue) - (Head Circumference: 21 – 23.25 inches)	\$99.00
2000FA-*	Foam Abductor	\$105.00
2000HICP-*	Hip Inflection Comfort Pommel	\$175.00
2000HIW-*	Hip Inflection Wedge	\$105.00
2000IC15-*	Quick-Change Incontinent Cover for 15 Inch Seat Depth Base	\$145.00
2000IC18-*	Quick-Change Incontinent Cover for 18 Inch Seat Depth Base	\$145.00
2000ADD	Adductor Strap	\$75.00
2000AB	Attachment Bracket Kit (for Vehicles without LATCH) – See Important Note Below****	\$95.00

<p>* The "Churchill Standard" includes the following features:</p> <ul style="list-style-type: none"> • Base • Lower Anchor Attachment Straps • Torso Support System • Anti-Slouching Leg Strap Assembly • Foam Padding • Fabric Cover • Positioning Vest (Small, Medium, or Large) • Tether Extension Assembly • Belt Positioning Guides <p>** One Vest is included in price of "Standard Churchill with Vest." Retail Price of additional vests is \$375.00 each.</p>	<p>***The "Churchill with Harness" includes the following features:</p> <ul style="list-style-type: none"> • Base • Lower Anchor Attachment Straps • Positioning Harness • Foam Padding • Fabric Cover • Tether Extension Assembly • Belt Positioning Guides
--	--

**** The Attachment Bracket Kit is designed for use in vehicles that do not have LATCH. The Churchill must be installed with LATCH or with the Attachment Bracket kit. All Vehicles 2003 and newer are equipped with LATCH in at least two rear seating positions, so there should be no installation issues for these vehicles and the Attachment Bracket will not be necessary. Many 2000 - 2002 models also have LATCH (See your vehicle owner's manual to verify if LATCH is present).

For vehicles without LATCH, the attachment bracket can only be used on vehicle seats that you can slide a ruler, or similar object, from the back side of the vehicle seat through the crack in the vehicle seat (where the seat bottom and seat back come together) through to the front side of the seat. This can usually be done in mini-vans and SUVs. Generally, the Attachment Bracket cannot be installed in the back seat of sedans. Installation of the Attachment Bracket Kit does require that a hole be drilled through the vehicle's floor board.

Shipping charges – The cost of shipping "Best Way" will be added to Purchase Orders.

Thanks for Your Business!