

Hope Car Bed

by Merritt Manufacturing



PRONE POSITION

Attached sleeping bag with adjustable cummerbund



RIGHT SIDE POSITION

Attached sleeping bag with adjustable cummerbund



SUPINE POSITION

3-point harness with adjustable cummerbund (for chest protection)




PRONE POSITION

3-point harness with adjustable cummerbund (for back protection)

SUPINE POSITION

Attached sleeping bag with adjustable cummerbund

Unique travel solutions for children with special health care needs requiring a larger car bed



FEATURES

- Easy to install securely in vehicle
- EPS Foam
- Inner and outer plastic shell
- Car bed shell 10 year life - all parts replaceable
- Installation training video available
- Easy to clean cushions/covers
- Rebound control

SPECIFICATIONS

- Child weight range 4.5 - 35 lbs.
- Child maximum height 29" - longer if knees are permitted to bend
- Inside bed dimensions: 14.5" X 29.5"
- 5 Restraint positions
- Small and large size sleeping bags plus side load bag available
- 3-Point restraint option
- Car bed weight 22.2 lbs.