

The Jefferson Car Seat for Omphalocele Patients



Now available at Merritt Manufacturing Inc.

The Jefferson Design Features

- Rear Facing 7.5 -40 lbs
- Support cushion for the smaller child (two sizes available)
- Durable medical grade naugahyde cleanable cover
- Safeguard latch connectors



Rear Facing Seat 7.5-40 lbs

Design Features

- Rebound Control Bar – prevents child from contacting vehicle seat on rebound
- E-Z Leveling feature– provides secure adjustment to the medically desired angle up to 50 degrees
- Expanded leg room
- Improved access to harness adjustment strap



E-Z Leveling feature – Eliminates the need for rolled blankets and/or pool noodles

Design Features



- Yoke Harness design – Easy to use
- Allows trach to be safely used
- Two buckle positions – prevents buckle/harness contact with the omphalocele
- Soft foam buckle pad to protect child from buckle

Yoke Harness Design