



2022 Wallenberg Retail Price List

The Wallenberg-Pediatric Positioning
Vehicle Restraint System

Part Number	Description	Price
6000WS-N	<p>The Wallenberg Standard* a pediatric vehicle restraint system for children with hip or arm casts</p> <p>Weight range: 5 to 80 lbs.;</p> <p>Rear Facing Weight Range: 5 to 40 lbs.;</p> <p>Forward Facing Weight Range: 25 to 80 lbs. (We strongly recommend rear facing to the top rear facing weight and height range of the Wallenberg).</p>	\$2,710.00

* The "Wallenberg Standard" includes the following features:

- Absolutely no sides, which enables this seat to completely accommodate hip casts
- Also designed for arm casts
- Most webbing is removable for easy disinfecting
- Cover and foam are disposable
- A "Hammock" designed specifically for children in Hip Casts
- Stainless steel and plastic construction make it extremely durable and easy to clean

Shipping charges – The cost of shipping "Best Way" will be added to Purchase Orders.