



Fall 2020

nova

Bathroom Safety





LESS THAN A DOLLAR A DAY

We use the bathroom every day ... 365 days of each year, yet this room is often the most dangerous place in your home. Being proactive with effective product solutions is the key to fall prevention.

For less than a dollar a day, you can equip your bathroom with products to keep you safe. Make your bathroom happy and safe every day.

Do You Know the **BATHROOM HOT SPOTS?**

GET SOME LEVERAGE



HOT SPOT #1

If you have trouble getting on and off a chair, then you are also having trouble getting on and off a toilet. An elevated toilet seat and hand rails can make all the difference.

GETTING IN & OUT OF THE TUB



HOT SPOT #2

Getting in and out of the tub or shower can be tricky, and holding onto the sink or towel rack isn't safe. Using grab bars is a much better solution.

BATHING INDEPENDENTLY



HOT SPOT #3

Our tub or shower is our sanctuary, but without the proper safety equipment, this is another bathroom hot spot. Bath seats, hand-held shower heads and bath mats are a must!

Bath Stool
9006-R



- Great for travel and easy to store
- Tool-free assembly in minutes
- Adjustable height aluminum frame
- Latex-free rubber tips

| | |
|---|---|
| • Seat Height Adjustment: 15.5" - 20.5" | • Leg Footprint: High - 15"w x 14.5"d |
| • Seat Dimensions: 12.5"w x 12.75"d | • Weight Capacity: 300 lb |

Bath Stool
9005



- Great for use in small areas
- Adjustable height aluminum frame
- Latex-free rubber tips

| | |
|--|--|
| • Seat Height Adjustment: 13.75" - 20" | • Leg Footprint: High - 15.5"w x 14.75"d |
| • Seat Dimensions: 12.75"w x 12.75"d | • Weight Capacity: 250 lb |

Tub Shower Board
9200

- Width adjusts to attach securely to most bathtubs
- Built-in soap holder
- Rubber covered clamps resist slipping and protect surfaces from scratching
- Safety handle can be attached to either side for support



| | |
|--|----------------------------------|
| • Seat Dimensions: 31"w x 14.5"d | • Weight Capacity: 250 lb |
| • Width Adjustment: 17.75" - 28.5" | |

Shower Chair
9401 with U-Shaped Cutout



- U-shaped design for cleaning purposes
- Arms provide support for sitting and rising
- Adjustable height aluminum frame
- Latex-free rubber tips

| | |
|--|---|
| • Seat Height Adjustment: 18" - 20.75" | • Leg Footprint: High - 21.25"w x 17.25"d |
| • Seat Dimensions: 16"w x 15.75"d | • Weight Capacity: 250 lb |

Foldable Shower Chair
9400 with U-Shaped Cutout



- Great for travel/easy to store
- Folds and unfolds quickly
- U-shaped design for cleaning purposes
- Arms provide support for sitting and rising
- Adjustable height aluminum frame
- Latex-free rubber tips

| | |
|--|--|
| • Seat Height Adjustment: 18" - 20.75" | • Leg Footprint: High - 21"w x 19.25"d |
| • Seat Dimensions: 16"w x 16"d | • Weight Capacity: 250 lb |

Wall Mounted Shower Seat
9404

- Mounts directly to shower wall
- Saves space by folding up when not in use
- Includes wall mount anchors
- Adjustable height aluminum frame
- Latex-free rubber tips



| | |
|---|---|
| • Seat Height Adjustment: 19.75" - 22.25" | • Overall Footprint: High - 15.75"w x 13.75"d |
| • Seat Dimensions: 16"w x 16"d | • Weight Capacity: 250 lb |

Deluxe Transfer Bench**9070**

- Built-in soap holder and hand held shower holder
- Converts easily for left or right hand transfers
- Adjustable height aluminum frame
- Latex-free rubber tips & suction tips



- | | |
|---|---|
| • Seat Height Adjustment: 18.75" - 23.25" | • Leg Footprint: High - 33"w x 17.5"d |
| • Seat Dimensions: 31.5"w x 14.5"d | • Weight Capacity: 300 lb |

Standard Transfer Bench**9071**

- Converts easily for left or right hand transfers
- Adjustable height aluminum frame
- Latex-free rubber tips & suction tips



- | | |
|---|--|
| • Seat Height Adjustment: 17.25" - 21.25" | • Leg Footprint: High - 30"w x 20.75"d |
| • Seat Dimensions: 27"w x 16.25"d | • Weight Capacity: 300 lb |

Heavy Duty Transfer Bench**9075**

- Converts easily for left or right hand transfers
- Adjustable height aluminum frame
- Latex-free rubber tips & suction tips



- | | |
|---|--|
| • Seat Height Adjustment: 16.75" - 21.75" | • Leg Footprint: High - 29.5"w x 20.75"d |
| • Seat Dimensions: 27"w - 16.25"d | • Weight Capacity: 500 lb |

Padded Transfer Bench**9080**

- Padded seat and back add comfort when bathing
- Converts easily for left or right hand transfers
- Adjustable height aluminum frame
- Latex-free rubber tips & suction tips



- | | |
|---|---|
| • Seat Height Adjustment: 17.75" - 21.75" | • Leg Footprint: High - 29.75" x 20.5"d |
| • Seat Dimensions: 27"w x 16"d | • Weight Capacity: 300 lb |

Portable Transfer Bench**9072-R**

- Great for travel and easy to store
- Tool-free assembly in minutes
- Converts easily for left or right hand transfers
- Latex-free rubber tips & suction cups



- | | |
|---|--|
| • Seat Height Adjustment: 19.25" - 23.25" | • Leg Footprint: High - 31.75"w x 22.5"d |
| • Seat Dimensions: 26.5"w x 16.25"d | • Weight Capacity: 300 lb |

Padded Transfer Bench/Commode**9073**

- Converts easily for left or right hand transfers
- Removable seat for easy cleaning
- Bucket and lid included
- Adjustable height aluminum frame
- Latex-free rubber tips & suction tips



- | | |
|---|--|
| • Seat Height Adjustment: 19" - 23.5" | • Leg Footprint: High - 36.5"w x 17.75"d |
| • Seat Dimensions: 26.75"w x 14.5"d | • Weight Capacity: 300 lb |

Bedside Commode

8450



- Can be used beside the bed and over the toilet as a raised toilet seat with arms
- Adjustable height steel frame
- Bucket, lid, and splash guard included
- Toilet seat opening: 8.5" w x 9.5" d

- | | |
|---|---|
| • Seat Height: 17.5" - 23.5" | • Seat Dimensions: 13.25" w x 14.25" d |
| • Overall Dimensions: 24" w x 19.75" d | • Weight: 11.5 lb |
| • Width Between Arms: 17.75" | • Weight Capacity: 350 lb |

Folding Commode

8700-R



- Can be used beside the bed and over the toilet as a raised toilet seat with arms
- Ideal for travel: folds and unfolds in seconds
- Adjustable height steel frame
- Bucket, lid, and splash guard included
- Toilet seat opening: 8" w x 10" d

- | | |
|--|---|
| • Seat Height: 16.25" - 22.25" | • Seat Dimensions: 13.75" w x 16.25" d |
| • Overall Dimensions: 22.5" w x 21.5" d | • Weight: 10 lb |
| • Width Between Arms: 18.5" | • Weight Capacity: 300 lb |

Drop-Arm Commode

8900W



- Arms drop for easy patient transfer
- Bucket, lid, and splash guard included
- Available in Gray and Chrome
- Adjustable height steel frame
- Toilet seat opening: 7.75" w x 9.75" d

- | | |
|--|---|
| • Seat Height: 18.75" - 22.25" | • Seat Dimensions: 13.5" w x 16" d |
| • Overall Dimensions: 22.5" w x 20.75" - 22.75" d | • Weight: 19.6 lb |
| • Width Between Arms: 18" | • Weight Capacity: 300 lb |

Padded Drop-Arm Commode

8901W



- Arms drop for easy patient transfer
- Padded seat and backrest
- Bucket, lid, and splash guard included
- Adjustable height steel frame
- Toilet seat opening: 8.25" w x 11" d

- | | |
|---|---|
| • Seat Height: 19.25" - 23.75" | • Seat Dimensions: 16.5" w x 16.5" d |
| • Overall Dimensions: 22.5" w x 21.25" - 23" d | • Weight: 20 lb |
| • Width Between Arms: 17.25" | • Weight Capacity: 300 lb |

Heavy Duty Commode

8580



- Extra-large seat provides added support
- Bucket and lid included
- Adjustable height steel frame
- Toilet seat opening: 9" w x 11" d

- | | |
|---|--|
| • Seat Height: 19" - 23" | • Seat Dimensions: 22.75" w x 18.5" d |
| • Overall Dimensions: 34.75" w x 25" d | • Weight: 22.25 lb |
| • Width Between Arms: 26.75" | • Weight Capacity: 450 lb |

Heavy Duty Commode

8582



- Extra-large seat provides added support
- Bucket, lid, and splash guard included
- Adjustable height steel frame
- Toilet Seat Opening: 8.25" w x 10.25" d

- | | |
|---|---|
| • Seat Height: 17.5" - 24.25" | • Seat Dimensions: 23.25" w x 18.75" d |
| • Overall Dimensions: 27" w x 23.75" d | • Weight: 16.6 lb |
| • Width Between Arms: 23.75" | • Weight Capacity: 500 lb |

Heavy Duty Drop-Arm Commode

8583

- Arms drop for easy patient transfer
- Extra-large seat provides added support
- Bucket and lid included
- Adjustable height steel frame
- Toilet Seat Opening: 8.25" w x 10.25" d



- | | |
|-------------------------------------|--|
| • Seat Height: 19.75" - 23.75" | • Seat Dimensions: 23.25" w x 18.75" d |
| • Overall Dimensions: 30" w x 22" d | • Weight: 18.1 lb |
| • Width Between Arms: 25.75" | • Weight Capacity: 500 lb |

Drop-Arm Transport Chair Commode

8805

- Arms drop for easy patient transfer
- 4" swivel wheels with locking casters
- Padded seat fits over commode and converts into transport chair
- Bucket & lid included
- Removable, swing away footrests
- Bucket can be removed without lifting user or raising seat
- Toilet Seat Opening: 8.75" w x 10" d
- NOT FOR USE IN SHOWER



- | | |
|--|---------------------------------------|
| • Seat Height: 18.75" | • Seat Dimensions: 17.75" w x 17.5" d |
| • Overall Dimensions: 21" w x 22.25" d x 38" h | • Weight: 22 lb |
| • Width Between Arms: 17.5" | • Weight Capacity: 250 lb |

Shower Chair/Commode

8800

- Waterfall armrests and padded seat
- 4" swivel wheels with locking casters
- Bucket and lid included
- Toilet seat opening: 8.5" w x 10" d



- | | |
|--|---------------------------------------|
| • Seat Height: 22.25" | • Seat Dimensions: 15.5" w x 16.75" d |
| • Overall Dimensions: 22.75" w x 22.5" d x 34" h | • Weight: 13.25 lb |
| • Width Between Arms: 19.5" | • Weight Capacity: 250 lb |

Deluxe Shower Chair/Commode

8801

- Waterfall armrests and padded seat
- 4" swivel wheels with locking casters
- Bucket and lid included
- Removable, swing away footrests
- Toilet seat opening: 8.75" w x 10.5



- | | |
|--|----------------------------------|
| • Seat Height: 19.5" | • Seat Dimensions: 16" w x 16" d |
| • Overall Dimensions: 22.75" w x 24.5" d x 36" h | • Weight: 18.1 lb |
| • Width Between Arms: 18.5" | • Weight Capacity: 250 lb |

Commode Liner

PA-1080

- Medical grade bag with super absorbent pad that quickly transitions body fluids into a gel, eliminating odors, spills & splashes
- Provides convenience & safety for Patient, Caregiver & Family
- No more hand-cleaning of soiled pail
- No handling of feces or urine
- Convenient drawstring closure for easier cleanup
- Pad absorbs 17 oz. of body fluids



Padded Toilet Seat Risers

2629-R (2") & 2630-R (4")

- High density foam minimizes pressure on sensitive areas
- Attaches easily to toilet seat with four adjustable Velcro straps to reduce shifting
- Easy to clean vinyl surface
- Fits most toilet seats
- Dimensions:
 - 2629-R: 14" w x 15.75" d x 2" h
 - 2630-R: 14.5" w x 15.75" d x 4" h
- Weight Capacity: 250 lb



2629-R



2630-R