

DISTRIBUTOR/AGENT:	DATE OF ORDER:
	DATE OF DELIVERY:
	NAME OF THE ADDRESSEE/USER:



FEATURES

- Non-toxic epoxy varnished chromium-plated **frame**, in steel.
- Padded **pelvic support** adjustable in height and backwards-forwards.
- Padded **thoracic support** adjustable in height and removable.
- Double thickness cotton fabric **harness** adjustable in height thanks to belts and buckles in girth with velcro strips. It can be easily removed for washing.
- **Unpierceable wheels**. 17 cm, front pivoting wheels and 30 cm back fixed wheels.
- Semi-elastic **PVC bumpers** to protect from impacts.
- **Hand brakes**.

Tested by **DEKRA Stuttgart**
 and listed in the GERMAN HMV (Technical Aids Listing)
FDA registered - USA
 Marked **CE**

Please check off the interested products



MEASURES

- DYNAMICO-3 (medium)**
YELLOW frame
- DYNAMICO-4 (large)**
ANTHRACITE frame
- DYNAMICO-5 (for adults)**
ANTHRACITE frame

Attention!
 Read carefully the user and maintenance manual



TS*: thoracic supp.
 PS*: pelvic supp.

Size	Max. Cap. Load	A	B	C	D	E	F	G	H	Total Weight	Weight Bars
3	Kg. 85	cm 85-110	cm 97-124	cm 25	cm 30	cm 85	cm 94	cm 74	cm 120-155	Kg. 19,3	Kg. 6,2
4	Kg. 100	cm 93-120	cm 100-134	cm 30	cm 35	cm 110	cm 101	cm 78	cm 140-175	Kg. 20,4	Kg. 7,4
5	Kg. 115	cm 106-133	cm 116-146	cm 32 TS* / 41 cm PS*	cm 41	cm 120	cm 111	cm 84	cm 160-195	Kg. 22,7	Kg. 8,7

ADDITIONAL COMPONENTS

<p><input type="checkbox"/> 809 Arm support. Alternative to 813</p>	<p><input type="checkbox"/> 943 Arm straps for arm supports 809</p>	<p><input type="checkbox"/> 813 Adjustable handlebar. Alternative to 809</p>	<p><input type="checkbox"/> 810 Weight bars (sizes 2,3,4,5). Please see the table above to check the added weight for each size</p>	<p><input type="checkbox"/> 814E Leg divider</p>
-------------------------------------------------------------------------------------------------	-------------------------------------------------------------------------------------------------	----------------------------------------------------------------------------------------------------------	-----------------------------------------------------------------------------------------------------------------------------------------------------------------	---------------------------------------------------------