

GRILLO mini

front drive and posterior gait trainer



POSTERIOR
PT



FRONT
DRIVE
PTM



6 versions
2 models: front drive and posterior,
just one folding frame
precise and continuous adjustments
interchangeable components
2 ergonomic and independent
supports for trunk and pelvis

 design and quality
MADE IN ITALY since 1980



FRONT DRIVE

POSTERIOR

AVAILABLE VERSIONS	FRONT DRIVE			POSTERIOR		
	PT	PA	PTM	PT	PA	P
INTERCHANGEABLE COMPONENTS	INCLUDING: TRUNK SUPPORT PELVIC SUPPORT ERGONOMIC HARNESS HANDLEBAR	INCLUDING: PELVIC SUPPORT ERGONOMIC HARNESS ARM SUPPORTS	INCLUDING: TRANSPARENT TRAY WITH HANDLEBAR PELVIC SUPPORT ERGONOMIC HARNESS	INCLUDING: TRUNK SUPPORT PELVIC SUPPORT ERGONOMIC HARNESS REMOVABLE GRIP HANDLES	INCLUDING: PELVIC SUPPORT ERGONOMIC HARNESS ARM SUPPORTS	INCLUDING: PELVIC SUPPORT ERGONOMIC HARNESS REMOVABLE GRIP HANDLES
925 Trunk Support	●	○	○ only when transparent tray is removed	●	○	○
924 Pelvic Support	●	●	●	●	●	●
930 Ergonomic harness	●	●	●	●	●	●
926 Handlebar	●	○	NA	NA	NA	NA
809 Arm Supports	○	●	NA	○	●	○
937 Multi-adjustable Handles	NA	NA	○	NA	NA	NA
927 Removable grip handles	NA	NA	NA	●	○	●
894 Back Belt	NA	NA	○	NA	NA	NA
890 Leg Divider	○	○	○	○	○	○
923 Pushing handlebar for care givers	○	○	○	○	○	○
929 Pedals for parking brakes	○	○	○	○	○	○

● = INCLUDED ○ = ACCESSORY NA = NOT AVAILABLE

- ALUMINIUM FRAME, foldable for easy transportation.
- PADDED PELVIC and TRUNK SUPPORTS, independently adjustable in height, width and circumference; allow an easy access thanks to the back opening. They hold the user just like two hands and this gives him/her a sense of security and stability allowing the user to concentrate only on walking without any fear.
- GRADUATED BARS as a reference on height adjustments.
- PELVIC HARNESS, ergonomic, comfortable, adjustable and easily removable

for washing. It supports the users in case of necessity without generating painful pressures.

- WHEELS, light and puncture-proof: 2 front swiveling, Ø15 cm and 2 rear solid, Ø 17,5 cm, for indoor and outside use.
- DIRECTIONAL LOCKS on front wheels.
- ANTIREVERSE device, ADJUSTABLE BRAKES, parking brake on back wheels.
- SEMI-ELASTIC BUFFERS on front wheels to protect against bumps

FRONT DRIVE	POSTERIOR	USER DIMENSIONS:
<p>Inner circumference pelvic and trunk support from cm 53 to cm 79</p> <p>Weight kg. 13,5 MAX. LOAD kg. 35</p>	<p>Inner circumference pelvic and trunk support from cm 53 to cm 79</p> <p>Weight kg. 13,5 MAX. LOAD kg. 35</p>	<p>USER HEIGHT (C): min cm 75 max cm 100</p> <p>IMPORTANT: Make sure that measures A and B indicated in the figure do not exceed the suggested range indicated for the pelvic and trunk support</p>



COMPANY WITH QUALITY MANAGEMENT SYSTEM CERTIFIED BY = ISO 9001 =