

DISTRIBUTOR/AGENT:	DATE OF ORDER:
	DATE OF DELIVERY:
	NAME OF THE ADDRESSEE/USER:



FEATURES

- **Ergonomic and modular tricycle** that helps the user to keep a proper posture and makes the ride easier thanks to its design and components.
- **Lowered frame** for easy ascent and descent.
- **Front positioned bottom bracket** for easier pedaling.
- **Optimized pinion** to facilitate the pedaling in safety.
- **Adjustable handlebar** in height, backward-forward and reversible.
- **Wide and ergonomic handles**, to increase comfort.
- **Bell.**
- **Back lamp.**
- **Cranks adjustable** in lenght.
- **Front brake** that can be operated by the lever either on the right or on the left side of the handlebar. The distance between the brake lever and the handlebar is adjustable by screw.
- **Adjustable ergonomic saddle** in height, tilt and backward-forward
- **Rear storage basket.**
- **12" solid wheels** run on ball bearings. A mechanism on the rear wheels allows easy adjustment of the tension of the chain.



**Attention!**  
Read carefully the user and maintenance manual


Please check off the interested products:

MEASURES

FRAME COLOUR:


- Blue
- Red
- Pink

View from above




Storage basket  
Internal measures

<b>L</b>	41 cm	In 16
<b>M</b>	18 cm	In 7
<b>N</b>	6,6 cm	In 2,6
<b>Max. Load</b>	15 kg	lb 33



C1= C2-2cm



C2= leg length

	Max. load	Aid weight	A	B	C1*	D	E	F	Use height(H)
<b>size mini</b>	Kg. 30 lb. 66,1	Kg. 15,5 lb. 34,1	cm 45 in. 17,7	cm 62 in. 24,4	cm 39-53 in. 15,3-20,8	cm 12 in. 4,7	cm 103 in. 40,5	cm 36 in. 14,2	cm 90-110 in. 35,4-43,3

\*With standard saddle. You can further adjust the height of the saddle (3cm ±) acting on the clamp.

CHOOSE THE CONFIGURATION >

**Fixed Pinion**

**Free Wheel**

Please **CHOOSE THE POSITION OF THE BRAKE** on Standard Handlebar or on Ergonomic Handles 880.

Left  Right

In the absence of indications, the brake will be mounted on the right.

OPTIONAL COMPONENTS, ALSO CHOOSABLE SEPARATELY

**880 Ergonomic handles adjustable in inclination independently (1 pair).** Same size as the Handlebar.

Choose the position of the brake on Ergonomic Handles





**845 Forearm support**

Left  Right




**875 Wrist retainer with support and multi-adjustable glove independently.**

Left  Right




**842P Pelvic support with belt, adjustable in height.**

Also select the component 883





**842T Trunk support, adjustable in height, inclination and backward-forward with belt and joint.**

Also select the component 883





**883 Rod adjustable in height for Thoracic/Pelvic support (single rod).**


Please select even the type of support.

**877 Flat backrest**




**843 Abductive saddle with holder for supports**





**840 Pedal cages with straps (1 pair).**




**841 Pedal cage with length and width adjustable foot imprints and fastening straps, with support kit.**

Left  Right




**846 Ankle support adjustable in height, with support kit. Only with component 841.**

Left  Right



**923 Additional assistant handle.**




**963 Steering assistant handle adjustable in height.**



**881 Steering limiter.**



**882 Steering limiter with springs to return the steering to the straight position.**




**879 Front basket.**

