

Standy



The Standy is available in two sizes. It's non-slip footboard makes independent transfer secure. The anatomical work table is height adjustable and removable.



The Standy's stabilization is a result of consolidated experience, constant research on materials, and the attempt to optimize posture. The adjustability of this unit's chest support, knee pads, and heel rests will fit most users. The non-slip footboard makes transferring from sitting to standing smoother for the user.



Accessories:



865 Headrest



902 Front Footrest



857 Side Supports



861 Transport Wheels

Measurements in Inches

Standy	Size 3	Size 4
Max Capacity	198 lbs	264.5 lbs
User Height	43 - 62.5	53 - 70.5
Unit Weight	64 lbs	84 lbs
A Chest pad width	12	14.5
B Knee pad to base	8 - 20	11 - 23
C Pelvic support to base	25 - 36.5	28.5 - 44
D Chest support to base	33 - 47	37 - 51
E Back wheel to front wheel	30	33
F Total depth	35	40
G Left wheel to right wheel	27	29.5
H Total width	29	31
I Chest pad to back supt width	5 - 11	8 - 12.5
L Table height	4 - 9	5 - 10

