

**Meywalk Mobility Device
Order Form 2016**

Account Information

Company Name: _____
 Address: _____
 City, State Zip: _____
 Phone: _____
 Fax: _____
 Email: _____
 Purchase Order: _____

Mark For: _____
 Ship To: _____
 Address: _____
 City, State Zip: _____
 Phone: _____

Order Information

To insure proper Meywalk sizing, please include required measurements below

Miniwalk I

Seat Height: 10 1/2" - 19"
 Trunk Support Height: 20" - 29"
 Ring Size: 8 1/4" w x 8" d
 Weight Capacity: 90 lbs

106061 Miniwalk I 2,997.00 Ea
 Includes: Castor Swivel Locks. Color: Turquoise

Miniwalk II

Seat Height: 14 - 25"
 Trunk Support Height: 27-35"
 Ring Size: 8 1/4" w x 8" d
 Weight Capacity: 90 lbs

106062 Miniwalk II 2,997.00 Ea
 Includes: Castor Swivel Locks. Color: Turquoise

Miniwalk I Accessories

- | | | | |
|--------|---|--------|----|
| 106000 | <input type="checkbox"/> Nylon Hygiene Ring Cover, 70cm | 85.00 | Ea |
| 106070 | <input type="checkbox"/> Antitip support | 270.00 | Pr |
| 106075 | <input type="checkbox"/> Leg Separation Panel | 540.00 | Ea |
| 106080 | <input type="checkbox"/> Legguides | 342.00 | Pr |
| 106085 | <input type="checkbox"/> Hip Pads | 570.00 | Pr |
| 106089 | <input type="checkbox"/> No-Roll Back Attachment | 330.00 | Pr |
| 106096 | <input type="checkbox"/> Drag brakes | 270.00 | Pr |

Meywalk 2000 Small **New**

Seat Height: 19"-26"
 Trunk Support Height: 30 3/4"-38 1/2"
 Ring Size: 10"w x 10"d
 Weight Capacity: 150lbs

104500 Meywalk 2000 Small 3,795.00 Ea
 Color: Green

Meywalk 2000 Medium

Seat Height: 21"-29"
 Trunk Support Height: 33"-40"
 Ring Size: 12 3/4"w x 11 1/2"d
 Weight Capacity: 220lbs

104175 Meywalk 2000 Medium 3,795.00 Ea
 Color: Blue

Meywalk 2000 Large

Seat Height: 27" -37"
 Trunk Support Height: 41" - 51"
 Ring Size: 14 1/2" w x 13 3/4" d
 Weight Capacity: 220lbs

104173 Meywalk 2000 Large 3,795.00 Ea
 Color: Blue

Meywalk 2000 Accessories

- | | | | |
|--------|--|--------|----|
| 104000 | <input type="checkbox"/> Nylon Hygiene Ring Cover, 105cm | 85.00 | Ea |
| 104001 | <input type="checkbox"/> Nylon Hygiene Ring Cover, 120cm | 85.00 | Ea |
| 104002 | <input type="checkbox"/> Nylon Hygiene Ring Cover, 90cm | 85.00 | Ea |
| 104023 | <input type="checkbox"/> No-Roll Back Attachment | 330.00 | Pr |
| 104044 | <input type="checkbox"/> Leg separation plate | 540.00 | Ea |
| 104110 | <input type="checkbox"/> Height reducing fittings, -3" | 138.00 | Pr |
| 104166 | <input type="checkbox"/> Drag Brakes | 270.00 | Pr |
| 104167 | <input type="checkbox"/> Antitip support | 270.00 | Pr |
| 104168 | <input type="checkbox"/> Legguides | 342.00 | Pr |
| 104215 | <input type="checkbox"/> Hip Pads | 570.00 | Pr |
| 104997 | <input type="checkbox"/> Front Castor Swivel Locks | 120.00 | Pr |
| 10250 | <input type="checkbox"/> Comfort Seat for S, M and L | 225.00 | Ea |

Required Measurements:

A: Ring Height _____
 Measure user from floor to axilla to indicate ring height.

B: Inseam Height: _____
 Measure user from floor to inseam, with user standing. Indicates seat height.

C: Trunk Width: _____
 Measure user across trunk, at axilla height, side to side in a straight line, to indicate proper ring size.

