

**Account Information**

Company Name: \_\_\_\_\_  
 Address: \_\_\_\_\_  
 City, State Zip: \_\_\_\_\_  
 Phone: \_\_\_\_\_  
 Fax: \_\_\_\_\_  
 Email: \_\_\_\_\_  
 Purchase Order: \_\_\_\_\_

Mark For: \_\_\_\_\_  
 Ship To: \_\_\_\_\_  
 Address: \_\_\_\_\_  
 City, State Zip: \_\_\_\_\_  
 Phone: \_\_\_\_\_

**Order Information**

**To insure proper Meywalk sizing, please include required measurements below**

**Miniwalk I**

Seat Height: 10 1/2" - 19"  
 Trunk Support Height: 20" - 29"  
 Ring Size: 8 1/4" w x 8" d  
 Weight Capacity: 90 lbs

106061  Miniwalk I 3,497.00 Ea  
 Includes: Castor Swivel Locks. Color: Turquoise

**Miniwalk II**

Seat Height: 14 - 25"  
 Trunk Support Height: 27-35"  
 Ring Size: 8 1/4" w x 8" d  
 Weight Capacity: 90 lbs

106062  Miniwalk II 3,497.00 Ea  
 Includes: Castor Swivel Locks. Color: Turquoise

**Miniwalk Accessories**

106000	<input type="checkbox"/> Nylon Hygiene Ring Cover, 70cm	85.00	Ea
106070	<input type="checkbox"/> Antitip support	295.00	Pr
106075	<input type="checkbox"/> Leg Separation Panel	600.00	Ea
106080	<input type="checkbox"/> Legguides	375.00	Pr
106085	<input type="checkbox"/> Hip Pads	590.00	Pr
106089	<input type="checkbox"/> No-Roll Back Attachment	375.00	Pr
106096	<input type="checkbox"/> Drag brakes	295.00	Pr
11761	<input type="checkbox"/> Miniwalk Push/Drag Bar	375.00	Ea

**Meywalk 2000 Small \*\*New\*\***

Seat Height: 19"-26"  
 Trunk Support Height: 30 3/4"-38 1/2"  
 Ring Size: 10"w x 10"d  
 Weight Capacity: 150lbs

104500  Meywalk 2000 Small 4,345.00 Ea  
 Color: Green

**Meywalk 2000 Medium**

Seat Height: 21"-29"  
 Trunk Support Height: 33"-40"  
 Ring Size: 12 3/4"w x 11 1/2"d  
 Weight Capacity: 220lbs

104175  Meywalk 2000 Medium 4,345.00 Ea  
 Color: Blue

**Meywalk 2000 Large**

Seat Height: 27" -37"  
 Trunk Support Height: 41" - 51"  
 Ring Size: 14 1/2" w x 13 3/4" d  
 Weight Capacity: 220lbs

104173  Meywalk 2000 Large 4,345.00 Ea  
 Color: Blue

**Meywalk 2000 Accessories**

104000	<input type="checkbox"/> Nylon Hygiene Ring Cover, 105cm	85.00	Ea
104001	<input type="checkbox"/> Nylon Hygiene Ring Cover, 120cm	85.00	Ea
104002	<input type="checkbox"/> Nylon Hygiene Ring Cover, 90cm	85.00	Ea
104023	<input type="checkbox"/> No-Roll Back Attachment	375.00	Pr
104044	<input type="checkbox"/> Leg separation plate	600.00	Ea
104110	<input type="checkbox"/> Height reducing fittings, -3"	175.00	Pr
104166	<input type="checkbox"/> Drag Brakes	295.00	Pr
104167	<input type="checkbox"/> Antitip support	295.00	Pr
104168	<input type="checkbox"/> Legguides	450.00	Pr
104215	<input type="checkbox"/> Hip Pads	590.00	Pr
104997	<input type="checkbox"/> Front Castor Swivel Locks	145.00	Pr

**Required Measurements:**

**A: Ring Height** \_\_\_\_\_  
 Measure user from floor to axilla to indicate ring height.

**B: Inseam Height:** \_\_\_\_\_  
 Measure user from floor to inseam, with user standing. Indicates seat height.

**C: Trunk Width:** \_\_\_\_\_  
 Measure user across trunk, at axilla height, side to side in a straight line, to indicate proper ring size.

