

miniwalk



- Moves up and down during gait, allowing for a more normal heel strike, toe off and mid phase as the child moves forward
- Spring assist system to aide the child to “swing through” during gait
- Standard with castor swivel locks. Can be engaged for in line mobility
- Adjustable angle handlebar encourages more upright or forward trunk orientation
- **AN EARLY START ...** ages 1 1/2 to 7 years (refer to sizing table below)

Dynamic Gait Movement

The Miniwalk moves up and down with each step, allowing the user to come down into a heel strike and toe off. As the user steps forward, the spring lifts and assists with the swing through phase of gait

Wide range of accessories.....

- Anti tippers (req'd for no back and drag brakes)
- No back brakes
- Drag Brakes
- Midline Leg Separation Panel
- Lateral Hip Pads
- Castor Swivel Locks (standard)

How to measure

Child's axilla height

Child's seat height



Sizing	Mini 1	Mini 2
Axilla Height	20-29 "	27-35"
Seat Height	10.5-19"	14-25"
Max users weight	90 lb	90 lb

Finally... a Gait Trainer with Dynamic Movement to promote normal gait motion!

pacific rehab inc

Designs for function