

Account Information

Company Name: _____ Address: _____ City, State Zip: _____ Phone: _____ Fax: _____ Email: _____ Purchase Order: _____	Mark For: _____ Ship To: _____ Address: _____ State Zip: _____ Phone: _____
--	---

Order Information

Meywalk Mk4 Medium		Meywalk Mk4 Large	
Seat Height*:	23.6 - 33" (Inseam)	Seat Height*:	32" - 40" (Inseam)
Trunk Support Height*:	39.8 - 52.3" (Axilla)	Trunk Support Height*:	48" - 59" (Axilla)
Ring Size**:	41.5" - 47" Circumfrence (105cm)	Ring Size**:	47" - 51" Circumfrence (120cm)
Weight Capacity:	285lbs (130kg)	Weight Capacity:	285lbs
Frame Color:	Red	Frame Color:	Black


011200 Meywalk MK4 Medium 4,314.00 Ea 011250 Meywalk MK4 Large 4,314.00 Ea

*To insure correct sizing, please include the required measurements below.
 **Other ring sizes available: 70, 90, 105, 120, 130cm circumference.

Mk 4 Medium Accessories

011147 <input type="checkbox"/> Antitip support, for Medium	270 Pr	010192 <input type="checkbox"/> Drag Brakes	269 Pr
011150 <input type="checkbox"/> Leg Divider	810 Ea	010193 <input type="checkbox"/> Non Reverse Brake	330 Pr
011160 <input type="checkbox"/> Lateral Leg Guides	870 Pr	010189 <input type="checkbox"/> Castor Swivel Locks	145 Pr
013075 <input type="checkbox"/> Hip Pads	595 Pr		
011730 <input type="checkbox"/> Tray	575 Ea		

Measuring and sizing the Meywalk MK4



***Required User Information**

Axilla Height: _____

Inseam Height: _____

User Weight: _____

Notes: