

# PANTHERA BAMBINO



**THE PANTHERA BAMBINO** is an easy rolling, everyday wheelchair for children. The Bambino is easy to drive, easy to adapt, and easy to lift. It's quite simply an active chair for children.

Ideal for children ranging in ages 4 to 12 years old who need an easy-to-handle everyday chair.

Now also available with abduction front (available only in red frame). Standard chassis, available in 4 widths and two colours (red and black) with custom colour options.

## STANDARD FEATURES

- Chrome molybdenum chassis
- Footrest adjustable in height, angle & depth
- 5" front caster
- Rear wheels
- Balancing adjustable axle
- Carbon fibre rear axle
- Dual adjustable anti-tippers
- Light weight frame: starting at 10.758lbs (4.89kg)

## ACCESSORIES (Sold Separately)



Adjustable Push Handle



Side Guards with Arm Rest



Push Bar



Cushion



Free-wheel

## PANTHERA BAMBINO SPECIFICATIONS

<b>Seat width</b>	9.45" (24 cm)	10.63" (27 cm)	11.8" (30 cm)	13" (33 cm)
<b>Backrest height</b>	6.7 - 9.8" (17 - 25 cm) , 10.6 - 13.8" (27 - 35 cm)			
<b>Seat height (rear)</b>	15" (39 cm)			
<b>Seat height (front)</b>	16.1" (41 cm)			
<b>Seat angle</b>	9-11.8" (23 -30 cm) adjustable			
<b>Seat depth</b>	4°			
<b>Total length</b>	26" - 28.1" ( 66 - 71.5 cm)			
<b>Total height</b>	22.4 - 25.6" (57-65 cm)			
<b>Total width</b>	20.5" (52 cm)	21.65" (55 cm)	22.8" (58 cm)	24" (61 cm)
<b>Max User weight</b>	110 lb (50 kg)			
<b>Total weight</b>	17.86 lb (8.12 kg)	18.36 lb (8.35 kg)	18.87 lb (8.58 kg)	19.37 lb (8.80 kg)
<b>Transfer weight</b>	10.77 lb (4.89 kg)	11.31 lb (5.14 kg)	11.81 lb (5.37 kg)	12.31 lb (5.60 kg)

**panthera**

Distributed exclusively  
in North America by

