

Adjusted by:

E-mail:

Adjusted for:

Tel:

Date:

Total weight (g): **0**

The weight changes based on the choices made below.

**S3 short low** has a adjustable rear axle in carbon fibre. To find a balance that suits you the rear axle can be adjusted into different positions. Note desired measure at "Settings". S3 short has a seat height of 47 cm in the front and 43 at the rear, without cushion. To adjust the seat height you can chose other rearwheel sizes or a raising Kit with suitable castor fork from the List of options. **S3 short low** is 50 mm shorter and 25 mm lower than the S3.

Max userweight: 100 kg for seat width 33-39

Save the guide to the "desktop" and open it with **Adobe Reader**, for best function Design your wheelchair by clicking in the circles in this guide. Blue, underlined text is clickable to **show info** or image. **BB** = seatwidth, **HH**= backrest height

**Original design**

**Alternative options**

**Other options**

[info](#)

**1) Chassis-colour**

- Black** - Anodic Black
- Clear**

Optional colour - A higher price and longer delivery time for optional colour  
[Optional colour](#)



**2) Seat width (BB)** [article number](#)

- 33 cm
- 36 cm
- 39 cm

Art.no include: chassis, seat, calf strap, backrest lock

Clear

**3) Seat depth**

- Clear
- Seat depth 35 cm**  
Std for all seat width

- Adjustable seat depth 27-33 cm**  
For seat width 33-39 cm



**4) Backrest**

- Standard**  
523BBHH VER2
- Clear

**Height (HH)**

- 20 cm
- 25 cm
- 30 cm
- 35 cm
- 40 cm
- 45 cm
- Clear

**Tapered** [Info](#)

- 525BBHH VER2
- Curved back** With backward curving the anti-tip device is complemented with "Knob anti-tip device S3 low profile", art.no. 4353121

**Other backrest**

[List of options](#) art. \_\_\_\_\_



**5) Rear axle**

- Carbon standard**  
33023BB
- Clear

- Carbon +15 mm**  
33023BC\*

\* seat width BB +1

Use this rear axle when choosing arm support at Accessories and/or increased seat height

**Other** rear axle or seat height

[List of options](#) art. \_\_\_\_\_



**6) Castor and fork**

- 120 mm [art.no.](#)
- Clear

- 90 mm [art.no.](#)
- 150 mm [art.no.](#)

**7) Footrest**

- Footrest Titanium**  
36001BB
- Footrest U-shape small heel bar**  
3600134 (width 33)
- Clear

**Titanium, extended, reinforced**

- 36111 BB [needed for Freewheel](#)
- Foldable footplates**  
[art.no.](#)
- Heelstrap for footplates**  
[art.no.](#)
- Plastic plate**  
36010BB

**Other footrest**

[List of options](#) art. \_\_\_\_\_



**8) Brake**

- Brake S3**  
4630000
- Clear

- One hand brake right**  
[art.no.](#)

**Other brake**

[List of options](#) art. \_\_\_\_\_



**9) Rear wheel - without pushrim**

- 24" **standard**  
[art.no.](#)
- Clear

- 22" **standard** [art.no.](#) Rear wheel sizes
- 25" **Spinergy X** black spokes  
[art.no.](#)
- 26" **Spinergy X** black spokes  
[art.no.](#)

**Other rear wheel**

[List of options](#) art. \_\_\_\_\_



**Tire - choose tire**

- Schwalbe Right Run** 22",24",25",26"
- Clear

- Rugged** 24"
- Schwalbe **Marathon Plus** 24",25",26"

**PU** 24" (massive)



**10) Pushrim**

- 24" Titanium**  
190DD11 (x2 / wheelchair)
  - Clear
- DD = Wheel diameter

- Aluminum**  
1904404 x2 (04=24", 05=25", 06=26")
- Paragrip**  
19051DDBT x2
- Tetragrip**  
19052DDBT x2

**Maxgrip** 22", 24" (contains latex)  
1903060 x2 (60=24", 80=22")

## Accessories



### Side guard

- [Side guard S3](#) [Video](#)  
art.no.
- Clear



### Arm support

- [Arm support S3](#) [Video](#)  
art.no.



### Side guard Carbon

- [Side guard Carbon](#)  
art.no.



### Push handles

- [Height adjustable S3](#) [Video](#)  
4252060 VER2
- [Height adjustable S3 +10](#)  
4252063 VER2  
Max. userweight 100kg
- Clear



### Push handle Foldable

- [Foldable S3](#) [Video](#)  
4245000



### Push bar

- [Push bar S3](#)  
42400BB VER2  
Max. userweight 100kg



### Anti-tip

- [Anti-tip](#) std seat height  
4353000 [Video](#)

[Other seat height](#)



### Cushion Wedge cushion

- [2,5 cm](#)  
572BB35
- [5 cm](#)  
575BB35
- Clear
- [Wedge cushion](#)  
57710BB



### Spoke protectors

- [24" logo](#)  
4500000
- [24" no logo](#)  
4500011
- [22" logo](#)  
4500100
- [22" no logo](#)  
4500111
- Clear



### Other

- [Tube protection front](#)  
9900251 (2 pcs)
- [Tube protection rear](#)  
9900251 (1 pcs)
- [Crutch holder S3](#) \_\_\_\_\_ pcs  
9900301 (x1/crutch)
- [Backrest wedge](#) \_\_\_\_\_ pcs  
4710012

## Set-up adjustments



### Back rest angle

Note measure (A) for desired back rest angle.  
14-28 mm. Standard=20 mm

A= \_\_\_\_\_ mm



### Rear axle adjustment

Note measure (A) for desired balance.  
Standard=40 mm

A= \_\_\_\_\_ mm



### Footrest height [Footrest - Foot plates](#)

Note measure (A) and (B) for desired footrest.  
[Bambino foot plate with Feewheel - info](#)

A= \_\_\_\_\_ mm      B= \_\_\_\_\_ mm



### Transportation marking

- YES do mark the wheelchair with
- attachment points for transport in disabled service vehicle.

Art.no.9000152

NOTE! See [manual](#) for more information.



Extra