

# ELEVATION

*ultra*  
lightweight



- Dynamic Seating System designed to meet a multitude of everyday activities
- Instant and easy seat height adjustment becomes a natural part of your daily life
- 10" of on-the-fly seat-height adjustment
- 30° of on-the-fly back rest angle adjustment
- Ultralight aluminum construction
- Patented Elevation Stability System

**Pella**  
MOBILITY

# ELEVATION

## ELEVATION LETS YOU DO MORE

Designed to enable greater function, independence, and comfort, the Elevation wheelchair promotes activities of daily living. Its seat height and backrest recline angle adjustment capabilities allow the user to self-adjust their sitting posture to best suit their current activity. The Elevation is capable of up to 10 inches of on-the-fly seat height/dump adjustment and 30° of back recline.

Dynamic back angle adjustability

Dynamic seat height adjustability

Fold-down back

Ultralight aluminum construction

High performance Spinergy wheels option

Centre of gravity axle adjustment

Multiple caster options

Six exclusive Elevation frame colors to choose from. Custom colors also available. Contact us for details.



Caster Options



Adjustable Footrest



Wheel Options



Adjustable Backrest Height



Camber Adjustment

Seat Width	12"–19"
Seat Depth	14"–19"
Front Seat Height	18"–21"
Back Angle	85°–115° (on-the-fly adjustable)
Dynamic Seat Height	10" of range
Weight Capacity	250 lbs
Chair Weight	22 lbs/10 kg (incl wheels)

### Options

- Flexguard™ side guard
- Scissor style brakes
- Uni-tine caster forks
- Frog leg suspension forks
- Natural fit handrims