

# FUZE T50 JR

# effortless

tilt



- Short wheelbase provides excellent maneuverability
- ISO tested
- Transit-ready option
- Width adjustable
- Depth adjustable


**pdg**  
MOBILITY

WHEELCHAIRS FOR FUNCTION AND INDEPENDENCE

[WWW.PDGMOBILITY.COM](http://WWW.PDGMOBILITY.COM)

# FUZE T50 JR

JUNIOR TILT-IN-SPACE AT ITS BEST!



Nine standard frame colors to choose from and custom colors also available. Contact us for details.



Adjustable seat height



Adjustable axle position



Adjustable seat depth and back rest



Transit tie-downs



Adjustable seat width

**Seat width**, adjustable 12" – 18"  
**Seat depth**, adjustable 13" – 18"  
**Seat height**, adjustable 13" – 20"  
**Back angle**, adjustable 90° – 120°  
**Back height** 20" and 25"  
**Dynamic tilt** 50°

**Overall width**  
 12" & 16" rear wheels seat width + 7.5"  
 20", 22", 24" rear wheels seat width + 9"

**Overall Length**  
**Weight Capacity**  
**Chair Weight**  
**Transit Options Available**  
**Warranty**

24" (no front rigging)  
 180 lbs  
 31 lbs (without front rigging, wheels, anti-tippers, armrests)  
 ANSI/RESNA WC-19  
 PDG's exclusive TGIF warranty (Thank Goodness It Fits)  
 Lifetime guarantee on the frame.

PDG recommends all wheelchairs be fitted to the consumer under the guidance of a health care professional knowledgeable of the particular needs of the consumer and the limitations of the product.

**USA**  
 PDG Product Design Group  
 PMB # 173 - 1124 Fir Avenue  
 Blaine, WA  
 98230, USA

**CANADA**  
 PDG Product Design Group  
 Unit 103 - 318 East Kent Avenue South  
 Vancouver, BC  
 V5X 4N6

**T: +1 604 323 9220**  
 F: +1 604 323 9097  
 Toll Free (in N. America): 888-858-4422  
 Toll Free Fax (in N. America): 877-323-9093  
[www.pdgmobility.com](http://www.pdgmobility.com)