

PO#	Tag:	Dealer Name:	Dealer Contact:
Mailing Address:		Shipping Address:	
City:	State:	Zip Code:	City:
Tel:	Fax:		Tel:
<b>FUZE FRAME TYPE</b> 51491 Fuze T50JR: 50 Degree Manual Tilt Wheelchair E1233 \$2,750		Confirm Order Via E-mail: YES NO Send PDG e-News: YES NO E-Mail:	

FUZE FRAME OPTIONS			
51492	Standard Growing Base & Seat Frame <small>Max Weight Capacity 180 lbs.</small>		STD
51493	Transit Tie-Down System (TTDS) <small>(TTDS includes 4 factory installed transport brackets. TTDS has been tested according to ANSI/RENSA WC 19. See Owners Manual for complete details.)</small>	K0108	\$250

SEAT WIDTH - (Adjustable) (1" Diameter Tubing)			
51494	12" - 15" Seat Width		NC
	Set At: 12" 13" 14" 15"		
51498	15" - 18" Seat Width		NC
	Set At: 15" 16" 17" 18"		

SEAT DEPTH - (Adjustable)			
51502	13" Seat Depth		NC
51503	14" Seat Depth		NC
51504	15" Seat Depth		NC
51505	16" Seat Depth		NC
51506	17" Seat Depth		NC
51507	18" Seat Depth	E2203	NC

FRAME COLOR - (Seat Frame Grey)			
51508	Blue Shimmer		NC
51509	Tracer Red		NC
51510	Black Shimmer		NC
51796	Orange Copper		NC
51512	Misty Purple		NC
51619	Forest Green		NC
52860	Bubble Gum Pink		NC
52858	Orange Fusion		NC
52859	Berry Purple		NC
51513	Custom Color: _____ (Contact PDG for Specification)		\$300

ARM TYPE			
51514	No Armrests		STD
51515	Single Post Fixed Height Armrest - 12"		NC
51516	Single Post Adjustable Height Armrest - 10" - 14"	E0973	\$250
51517	Locking Flip-up Armrests <small>Not recommended for clients with excessive strength or tone. Not compatible with fold down back, curved armpad, back upholstery.</small>	E0973	\$250
	Set At: 8" 9" 10" 11" 12" 13" 14"		
	Dual Post Adjustable Height Armrest	E0973	\$295
51527	Option: 8" - 11"		
51528	Option: 10" - 14"		
51529	Option: 12" - 16"		

ARMPAD OPTIONS			
51530	Curved Armpad (pair) <small>Not compatible with Lap Tray, Locking Flip-up Armrests or Dual Post Armrests.</small>		STD
51531	Desk Length Armpad, 2" X 10" (pair)		NC
51532	Full Length Armpad, 2" X 14" (pair)		NC
	3" Wide Arm Pads <small>May interfere with wheels in some tilt configuration.</small>		\$105
51533	Full Length 3"x12" (pair)		
51534	Desk Length 3"x10" (pair)		
	Gel Armpad <small>May interfere with wheels in some tilt configuration.</small>		\$115
51535	Full Length 2"x12" (pair)		
51536	Full Length 3.5"x12" (pair)		
	Contour Arm Trough with Hand Pad <small>Not compatible with Lap Tray, Locking Flip-up Armrests or Dual Post Armrests. May interfere with wheels in some tilt configuration.</small>	E2209	\$125
53249	Left (4.5" X 18.5")		
53250	Right (4.5" X 18.5")		

SEAT FRAME OPTIONS			
51541	No Seat Pan		STD
51542	Solid Seat Pan	E2231	\$95

CASTER OPTIONS			
51544	5" x 1" Urethane		NC
51545	6" x 1-1/4" Urethane		NC
51546	6" x 1-1/4" Pneumatic		\$75
51547	6" x 1-3/4" Semi-Pneumatic	K0075	\$85
51548	8" x 1-1/4" Urethane		NC
51549	8" x 1-1/4" Pneumatic		\$75
51550	8" x 2" Urethane Caster	K0074	\$145
51551	8" x 2" Pneumatic Caster		\$145

PO#	Tag:	Dealer Name:	Dealer Contact:
-----	------	--------------	-----------------

REAR WHEEL OPTIONS		
51552	12" Composite Wheel	NC
51553	16" Composite Wheel	NC
51554	20" Composite Wheel	NC
51555	20" Spoke Wheel	\$75
51556	22" Composite Wheel	NC
51557	22" Spoke Wheel	\$75
51558	24" Composite Wheel	NC
51559	24" Spoke Wheel	\$75

WHEEL LOCKS, ANTI-TIPPERS		
51798	Foot Operated Wheel Locks <small>Not compatible with 1 arm drive.</small>	\$200
51567	Push To Lock Wheel Locks	STD
51568	Pull To Lock Wheel Locks	NC
51569	Attendant Wheel Lock <small>No Extension.</small>	\$99
51570	Wheel Lock Extensions <small>E0961</small>	\$60
51571	Rear Anti-Tippers (Pair) <small>Required safety option. E0971</small>	\$135

SEAT HEIGHT OPTIONS						
Check desired seat-to-floor height in chart below.						
<small>All dimensions are approximate and depend on rear wheel type and pre-tilt setting (shipped standard with 3° pre-tilt). Actual seat height can vary by as much as 1" but can be adjusted down by removing the pre-tilt collar/bumper.</small>						
Rear Wheel	5" Casters			Rear Wheel	6" Casters / 6" Fork	
12"	-	13"	14"	12"	14"	18" Inverted
16"	14"	15"	16"	16"	14"	15"
20"	16"	17"	-	20"	16"	17"
22"	17"	-	-	22"	17"	18"
				24"	18"	-
Rear Wheel	6" Casters / 8" Fork			Rear Wheel	8" Casters / 8" Fork	
12"	18" Inverted			12"	18" Inverted	
	19" Inverted				19" Inverted	
16"	-	15"	16"		20" Inverted	
20"	16"	17"	18"	16"	20" Inverted	
22"	17"	18"	19"	20"	16"	17"
24"	18"	19"	-	22"	17"	18"
				24"	18"	19"
						20"

BACK POST STYLE (7/8" Diameter Tubing)		
51572	Push Bar Back, Straight	STD
51573	Cane Handle Back, Straight <small>Not compatible with back upholstery.</small>	NC
51574	Push Bar Back, 8° Bend	\$75
51575	Cane Handle Back, 8° Bend <small>Not compatible with back upholstery.</small>	\$75

BACK POST HEIGHT (Seat Pan to Handle)		
51576	20"	STD
51577	25"	NC

BACK POST OPTIONS		
51578	Angle/Height Adjustable Push Bar <small>Not compatible with Cane Handle back.</small>	\$235
51579	Removable Cane Handle Push Bar <small>Not compatible with Push Bar back.</small>	\$75
	<b>Seat To Back Angle:</b>	STD
51580	Set At: 90° 96° 102° 108° 114° 120°	
51586	Fold Down Back	STD
52289	Fixed Back	\$75

REAR TIRE OPTIONS		
51560	Urethane Tires	STD
51561	Pneumatic Tires <small>Pneumatic tires may result in a finished seat height 1/2" higher than stated.</small>	\$60
51562	Pneumatic Tires With Flat-Free Inserts <small>Pneumatic tires may result in a finished seat height 1/2" higher than stated. Flat Free Inserts not compatible with Spoked Wheels.</small>	E2213 \$95

HEADREST OPTIONS		
51588	Headrest Mounting Plate <small>Not compatible with head rest extension, Angle/height adjustable pushbar, Cane handle back. Not required if mounting to a solid back.</small>	\$75

AXLE TYPE		
51563	Quick Release Rear Axles	STD
51564	Fixed Rear Axles	NC

TILT / RECLINE OPTIONS		
51620	Tilt Limiting Device - Limits Tilt At < 50°  Set Tilt Limit At: 20° 25° 30° 35° 35° 40° 45° 50°	STD
51595		
53143	30° Reclining Backrest <small>Not compatible with headrest extension, cane handle back, spring back suspension, fold down back or tension adjustable upholstery. Not compatible with 13" Seat Depth.</small>	\$830

HAND RIM OPTIONS		
51635	Omit Hand Rims	NC
51565	Standard Hand Rims	STD
51566	Plastic Coated Hand Rims	\$95

FRONT RIGGING		
51602	Omit Front Rigging	NC
51603	Sector Blocks Only (No Front Rigging Supplied)	NC
52290	External Front Rigging Mount	\$75
51604	Swing Away 70° Front Rigging	NC
51605	Swing Away 90° Front Rigging (3"-9") <small>Not compatible with 8" casters.</small>	STD
51606	Swing Away Elevating Legrests <small>Not available on 12" or 13" wide.</small>	E0990 \$355



# FUZE T50 JR

**USA Price List**

Effective January 31<sup>st</sup>, 2016 Page 3 of 4

PO#

Tag:

Dealer Name:

Dealer Contact:

## FOOT PLATES

51607	Composite Footplates Not compatible with 12" wide chair.		STD
51608	Adjustable Angle Footplates	K0040	\$175
51609	One Piece Flip-up Footplates		\$395
51839	Interlocking Angle Adjustable Footplate 14"-21" wide only		\$195

## EXTENSION TUBE

51613	Standard Length (13" - 17") Not compatible with 90°		STD
51614	Long Extension Tubes (15" - 19") Not compatible with 90°	K0043	\$75
51784	90° Extension Tubes (11" - 16") Not compatible with 60°, 70°, ELR.		NC

## FOOTREST/LEGREST ACCESSORIES

51610	Adjustable Length Heel Loops	E0951	\$50
51611	Body Point Calf Strap	K0038	\$60
52291	Swing Away 70° Calf Pad Not available on 12" or 13" wide.		\$135
52292	Miller's Dynamic footrest Mounts for Swing Away 70° Front Rigging		\$380
52293	Swing Away 70° Gel Padding		\$99
51806	Swing Away Elevating Legrest Gel Knee Pad		\$99

## POSITIONING STRAPS

51618	BodyPoint - Chest Strap	E0960	\$75
BodyPoint - PADDED - PEDIATRIC			
53261	Auto Style - HIP belt 1"	E0978	\$90
53262	4 point - HIP belt	E0978	\$110

## MISCELLANEOUS OPTIONS

51797	Clear lap tray, slide on. 15" - 18" wide only.		\$490
-------	--	--	-------

PO#	Tag:	Dealer Name:	Dealer Contact:
-----	------	--------------	-----------------

## OPTIONS FORM



MATRIX

MATRIX

CUSHION OPTIONS - MATRIX			
<b>MATRIX VI</b>	<input type="checkbox"/> 14" w x 15" d <input type="checkbox"/> 14" w x 17" d		\$550
<input type="checkbox"/> 53109	<input type="checkbox"/> 16" w x 16" d <input type="checkbox"/> 16" w x 18" d		
	<input type="checkbox"/> 18" w x 16" d <input type="checkbox"/> 18" w x 18" d <input type="checkbox"/> 18" w x 20" d		

BACK OPTIONS - MATRIX			
<b>MATRIX ELITE DEEP</b>	<input type="checkbox"/> 14" - 15" w x 16" h		\$925
<input type="checkbox"/> 52967	<input type="checkbox"/> 16" - 17" w x 16" h <input type="checkbox"/> 16" - 17" w x 18" h		
	<input type="checkbox"/> 18" - 19" w x 16" h <input type="checkbox"/> 18" - 19" w x 18" h		
<input type="checkbox"/> 52978	<b>MATRIX ELITE HEADREST ADAPTER</b>		NC
	Required if ordering Matrix Elite back and Headrest		

Please note with some configurations Matrix Back Hardware may interfere with fold down back.



VARILITE

VARILITE

CUSHION OPTIONS - VARILITE			
<b>VARILITE REFLEX</b>	<input type="checkbox"/> 16" w x 16" d <input type="checkbox"/> 16" w x 18" d		\$375
<input type="checkbox"/> 52980	<input type="checkbox"/> 18" w x 16" d <input type="checkbox"/> 18" w x 18" d <input type="checkbox"/> 18" w x 20" d		

BACK OPTIONS - VARILITE			
<b>VARILITE ICON TALL</b>	<input type="checkbox"/> 13" w - 15" w <input type="checkbox"/> 17" w - 19" w		\$950
<input type="checkbox"/> 52987	<input type="checkbox"/> 14" w - 16" w		
	<input type="checkbox"/> 15" w - 17" w		
	<input type="checkbox"/> 16" w - 18" w		




*A ROHO® for every need.*

ROHO

ROHO

CUSHION OPTIONS - ROHO			
<b>ROHO QUADTRO SELECT HIGH PROFILE</b>	<input type="checkbox"/> 16" w x 16" d <input type="checkbox"/> 16" w x 18" d <input type="checkbox"/> 16" w x 20" d		\$690
<input type="checkbox"/> 52993	<input type="checkbox"/> 18" w x 16" d <input type="checkbox"/> 18" w x 18" d <input type="checkbox"/> 18" w x 20" d		

BACK OPTIONS - ROHO			
<b>ROHO AGILITY MID CONTOUR BACK (WITH AIR ONLY)</b>	<input type="checkbox"/> 15" - 16" w x 16" h <input type="checkbox"/> 15" - 16" w x 20" h		\$840
<input type="checkbox"/> 53001	<input type="checkbox"/> 17" - 18" w x 16" h <input type="checkbox"/> 17" - 18" w x 20" h		



STEALTH

STEALTH

STEALTH HEADREST HARDWARE			
<input type="checkbox"/> 51589	<b>Multi-Axis Hardware - HMO475</b>	E1028	\$150
	Must have a headrest mounting plate OR solid back support for mounting		
<input type="checkbox"/> 51590	<b>Link - TWB480</b>	E1028	\$325
	Must have a headrest mounting plate OR solid back support for mounting		

STEALTH HEADREST PAD			
<input type="checkbox"/> 51593	<b>Stealth Headrest Pad Size</b>	E0955	\$200
	Set At: 10" <input type="checkbox"/> or 14" <input type="checkbox"/>		