

# STELLAR GL

## manual

45° of tilt



- Improved clearance for foot propulsion
- Ideal for hand propulsion with accessible rear wheels
- Knees remain low in tilt – wheel under desks, into vans and access other environments with ease
- Transit-ready option
- Short wheelbase provides excellent maneuverability
- Features adjustable seat height, seat depth and back angle
- Anterior and posterior tilt set-ups available



WHEELCHAIRS FOR FUNCTION AND INDEPENDENCE

[WWW.PDGMOBILITY.COM](http://WWW.PDGMOBILITY.COM)

