



# STELLAR GL

**Order Form - USA**

Prices in U.S. Dollars. January 1, 2020

## Customer and Order Information

Purchase Order *	Tag	v.0801
Dealer/Provider Name *		
Purchaser Name and Contact Information *		
Shipping Address *		
User Information (if available)	Weight _____ lbs	Height _____ ft/in <i>*required</i>

## Frame Type

<input checked="" type="checkbox"/> 1901000	<b>Stellar GL 45° Manual Tilt Wheelchair</b>	E1161	\$ 2,995
<input type="checkbox"/>	<b>Consider the lightweight Stellar GLT with aluminum frame.</b> See separate order form. _____		
<input type="checkbox"/>	<b>Consider the Stellar Impact for highly agitated users and the Stellar HD for bariatric users.</b> See separate order forms.		

## Frame Options

<input checked="" type="checkbox"/> 1902000	<b>Standard Base and Seat Frame</b> 250 lbs Weight Capacity.		
<input type="checkbox"/> 1902200	<b>Extended Base with Rear Floor Jacks</b>	K0108	\$ 900
	Lifts Rear Wheels off ground for increased stability. Limited ground clearance. Not compatible with Attendant Wheel Lock, HD Anti-Tipper, IV Pole, O2 Holder.		
<input type="checkbox"/> 1902300	<b>Transit Tie Down System (TTDS)</b>		
	TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RESNA WC 19. See Owners Manual for complete details. Not compatible with Attendant Wheel Locks.		

## Seat Width

Overall chair width is +9.5" with 12/16" Rear Wheels and +10.5" with 20/22/24" Rear Wheels.

<b>Standard Seat Widths</b>			
<input type="checkbox"/> 1903014	14" Seat Width	<input type="checkbox"/> 1903015	15" Seat Width
<input type="checkbox"/> 1903017	17" Seat Width	<input type="checkbox"/> 1903018	18" Seat Width
<input type="checkbox"/> 1903020	20" Seat Width	<input type="checkbox"/> 1903021	21" Seat Width
<input type="checkbox"/> 1903016	16" Seat Width	<input type="checkbox"/> 1903019	19" Seat Width
<input type="checkbox"/> 1903022	22" Seat Width		
		E2201	\$ 350

## Seat Depth

<b>Standard Seat Depths (Adjustable)</b>			
<input type="checkbox"/> 1904016	16" Seat Depth	<input type="checkbox"/> 1904017	17" Seat Depth
<input type="checkbox"/> 1904018	18" Seat Depth	<input type="checkbox"/> 1904019	19" Seat Depth
<input type="checkbox"/> 1904020	20" Seat Depth		
		E2203	\$ 500
<b>Custom Seat Depths</b>			
<input type="checkbox"/> 1904021	21" Seat Depth		
<input type="checkbox"/> 1904022	22" Seat Depth		
		E2204	\$ 750

## Frame Color

<b>PDG Premium Colors</b>			
<input type="checkbox"/> 1905000	Tracer Red	<input type="checkbox"/> 1905001	Gloss Black Shimmer
<input type="checkbox"/> 1905003	Fine Texture Black	<input type="checkbox"/> 1905002	Silver
<input type="checkbox"/> 1905006	Blue Shimmer	<input type="checkbox"/> 1905004	Orange Fusion
<input type="checkbox"/> 1905009	Geysers Green	<input type="checkbox"/> 1905005	Karl's Imagination
<input type="checkbox"/> 1905012	Bubble Gum Pink	<input type="checkbox"/> 1905007	Gloss White
<input type="checkbox"/> 1905200	<b>Custom Color</b> Contact PDG for specifications.	<input type="checkbox"/> 1905008	Dark Titanium
		<input type="checkbox"/> 1905010	Lollypop Blue
		<input type="checkbox"/> 1905011	Lollypop Red
		<input type="checkbox"/> 1905012	Berry Purple
		<input type="checkbox"/> 1905014	Daffodil Yellow
			\$ 300

### Armrest Type

<input type="checkbox"/> 1906000	<b>No Armrests</b>			
<b>Single Post Height Adjustable Armrests</b>			E0973	\$ 275
<input type="checkbox"/> 1906111	8"-11" Range	<input type="checkbox"/> 1906114	10"-14" Range	
<b>Dual Post Height Adjustable Armrests</b>			E0973	\$ 295
<input type="checkbox"/> 1906211	8"-11" Range	<input type="checkbox"/> 1906214	10"-14" Range	<input type="checkbox"/> 1906216
<b>Fixed Height Armrests (12")</b>			Full Length Armpad only.	
<input type="checkbox"/> 1906100	Single Post	<input type="checkbox"/> 1906200	Dual Post	
	<b>Armrest Position</b> Indicate preference for Armrest position. If left blank, the default setting will be used. In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.			
<input type="checkbox"/> 1906401	Rear	<input type="checkbox"/> 1906402	Mid	<input type="checkbox"/> 1906403
<b>Locking Flip-up</b>			E0973	\$ 295
Not recommended for users with excessive strength or tone. Not compatible with Fold Down Back, Clothing Guards, Back Upholstery. Requires Regular Full Length Armpads.				
<input type="checkbox"/> 1906310	Set at 10"	<input type="checkbox"/> 1906311	Set at 11"	<input type="checkbox"/> 1906312
<input type="checkbox"/> 1906313	Set at 13"	<input type="checkbox"/> 1906314	Set at 14"	

### Armpads & Armrest Options

<b>Regular Armpads</b>				
<input type="checkbox"/> 1907100	Desk Length - 2" x 10"	<input type="checkbox"/> 1907101	Full Length - 2" x 14"	
<b>Wide Armpads</b>			May interfere with wheels in some tilt configurations.	\$ 200
<input type="checkbox"/> 1907110	Desk Length - 3" x 10"	<input type="checkbox"/> 1907111	Full Length - 3" x 12"	
<b>Gel Armpads</b>			May interfere with wheels in some tilt configurations.	\$ 135
<input type="checkbox"/> 1907120	Full Length - 2" x 12"	<input type="checkbox"/> 1907121	Full Length - 3.5" x 12"	
<b>Contour Arm Trough with Hand Pad</b>			May interfere with wheels in some tilt configurations.	E2209
<input type="checkbox"/> 1907130	Right (4.5" x 18")	<input type="checkbox"/> 1907131	Left (4.5" x 18")	\$ 125
<input type="checkbox"/> 1907220	<b>Armrest Wideners for Dual Post Armrests</b> Limits tilt. Not compatible with 20", 22" or 24" Rear Wheels, Clothing Guards.			\$ 375
<b>Clothing Guards</b>			Not compatible with Flip-up Armrests, Armrest Wideners.	
<input type="checkbox"/> 1907310	For Single Post Armrests	Available for 10"-14" Height range only.		\$ 110
<input type="checkbox"/> 1907320	For Dual Post Armrests			\$ 165
<b>Armrest Extensions</b>			Not compatible with Flip-up or Single Post Armrests.	
<input type="checkbox"/> 1907400	3" Extension - Forward	<input type="checkbox"/> 1907401	3" Extension - Rearward	\$ 75

### Seat Pan

<input type="checkbox"/> 1908000	<b>No Seat Pan</b>	
<input type="checkbox"/> 1908001	<b>Solid Aluminum Seat Pan</b>	E2231 \$ 120

### Seat Height

**Don't forget** **Select desired Seat-to-Floor Height** STF Height is adjustable within the range of the selected Rear Wheel size.

RW	6" Casters					
12"	13"	14"	15"	16"		
16"			15"	16"	17"	18"
20"	13"	14"	15"	16"		
22"		14"	15"	16"	17"	
24"			15"	16"	17"	18"

RW	8" Casters					
16"			15"	16"	17"	18"
24"			15"	16"	17"	18"

Wheel size \_\_\_\_\_  
Caster size \_\_\_\_\_  
STF height \_\_\_\_\_

### Castors

<b>Urethane Castors</b>			8" Castors are compatible only with 16" and 24" Rear Wheels.	
<input type="checkbox"/> 1911061	6" x 1.25" Caster size	<input type="checkbox"/> 1911081	8" x 1.25" Caster size	
<input type="checkbox"/> 1911063	6" x 1.50" Caster size			\$ 85
<input type="checkbox"/> 1911084	8" x 2.00" Caster size			\$ 135



# STELLAR GL

Order Form - USA

Prices in U.S. Dollars. January 1, 2020

<b>Pneumatic Casters</b> 8" Casters are compatible only with 16" and 24" Rear Wheels.		
<input type="checkbox"/>	1911161 6" x 1.25" Caster size	E2214 & E2215 \$ 110
<input type="checkbox"/>	1911181 8" x 1.25" Caster size	E2214 & E2215 \$ 120
<input type="checkbox"/>	1911184 8" x 2.00" Caster size	E2214 & E2215 \$ 150
<input type="checkbox"/>	1911260 <b>Frog Legs 6" x 1.40" Soft Roll Casters</b>	\$ 95
<input type="checkbox"/>	1911400 <b>Frog Legs Suspension Forks</b> Not compatible with 12", 20", or 22" Rear Wheels. Requires 6" Soft Roll Caster.	E1015 \$ 395

## Rear Wheels

<b>Composite Wheels</b>						
<input type="checkbox"/>	1912012 12" Wheel size	<input type="checkbox"/>	1912016 16" Wheel size	<input type="checkbox"/>	1912020 20" Wheel size	
<input type="checkbox"/>	1912022 22" Wheel size	<input type="checkbox"/>	1912024 24" Wheel size			
<b>Spoke Wheels</b>					\$ 75	
<input type="checkbox"/>	1912120 20" Wheel size	<input type="checkbox"/>	1912122 22" Wheel size	<input type="checkbox"/>	1912124 24" Wheel size	
<input type="checkbox"/>	1912224 <b>24" Heavy Duty Spoke Wheels</b>					\$ 350
<b>Rear Wheel Position</b> Indicate Rear Wheel positioning preference. If left blank, the default setting will be used. In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.						
<input type="checkbox"/>	1912401 Rear	<input type="checkbox"/>	1912402 Mid	<input type="checkbox"/>	1912403 Front	

## Axle Type

<input type="checkbox"/>	1913000 <b>Quick Release Rear Axles</b>
<input type="checkbox"/>	1913010 <b>Fixed Rear Axles</b>

## Rear Tires

<input type="checkbox"/>	1914000 <b>Urethane Tires</b>	
<input type="checkbox"/>	1914010 <b>Pneumatic Tires</b> May result in a finished Seat Height 1/2" greater than stated.	E2211 & E2212 \$ 60
<input type="checkbox"/>	1914020 <b>Pneumatic Tires with Flat-Free Inserts</b> May result in a finished Seat Height 1/2" greater than stated. Not available with Spoke Wheels.	E2211 & E2213 \$ 95

## Hand Rims

Handrims are only available with 20", 22", and 24" Rear Wheels.

<input type="checkbox"/>	1915000 <b>No Handrims</b>	
<input type="checkbox"/>	1915010 <b>Standard Aluminum Handrims</b>	
<input type="checkbox"/>	1915020 <b>Plastic Coated Handrims</b>	\$ 135

## Wheel Locks

<b>User Wheel Locks</b>				
<input type="checkbox"/>	1916010 Push to Lock	<input type="checkbox"/>	1916020 Pull to Lock	
<input type="checkbox"/>	1916030 <b>Attendant Wheel Locks (Push to Lock)</b>			\$ 85
<input type="checkbox"/>	1916040 <b>User &amp; Attendant Wheel Locks</b> Available with 20", 22", or 24" Rear Wheels only.			\$ 125
<input type="checkbox"/>	1916100 <b>Wheel Lock Extensions</b> Not compatible with Attendant Wheel Lock.	E0961	\$ 60	
<input type="checkbox"/>	1916200 <b>Hill Holders</b> Requires Pneumatic Rear Tires with Flat-Free Inserts.	E0974	\$ 450	
<input type="checkbox"/>	1916300 <b>Drum Brakes</b> Available with 20" Rear Wheels only. <b>COMING SOON</b>	E2228	\$ TBD	

## Anti-Tippers

<input checked="" type="checkbox"/>	1917010 <b>Rear Anti-Tippers</b> Required Safety Option.	E0971	\$ 60
<input type="checkbox"/>	1917020 <b>Front Anti-Tippers</b>	E0971	\$ 795
<input type="checkbox"/>	1917030 <b>Heavy Duty Rear Anti-Tipper</b>	E0971	\$ 895

Provides additional rearward stability for users with excessive movement. Not compatible with Attendant Wheel Locks, IV Pole, O2 Holder.



# STELLAR GL

**Order Form - USA**

Prices in U.S. Dollars. January 1, 2020

## Back Post Style

- 1918010 **Push Bar Back**
- 1918020 **Cane Handle Back** Not compatible with Headrest Mounting Plate.
- Back Post Height** Please note: We have updated the descriptions/measurements below. The actual back post heights remain unchanged.
- 1918425 Short (22") *Previously 25"*  1918431 Tall (25") *Previously 31"* Measured from the Seat Pan to the Handle.

## Back Options

- Angle/Height Adjustable Push Bar Handle** Not compatible with Headrest Mounting Plate or Cane Handle Back. \$ 150
- 1919010 Button on Outside  1919011 Button on Inside
- Seat to Back Angle**
- 1919090 Set at 90°  1919098 Set at 98°  1919106 Set at 106°
- 1919114 Set at 114°  1919122 Set at 122°
- 1919210 **10" Backrest Extension** Requires the Angle/Height Adjustable Push Bar Handle. *E0966* \$ 150
- 1919212 **Headrest Mounting Plate** \$ 95  
Not compatible with Backrest Extension, Angle/Height Adjustable Push Bar, Cane Handle Back. Not required if mounting to a solid Back.
- 1919220 **Fold Down Back** Not compatible with Reclining Back.
- 1919230 **Spring Back Suspension** Not compatible with Reclining Back. \$ 105
- 1919240 **Seating Dynamics Dynamic Rocker Back Interface** *K0108* \$ 2,995  
Not compatible with Fold Down Back, Spring Back Suspension, 30° Recline, Cane Handle. Requires Adjustable Tension Upholstery or rigid Backrest. Finished Back Post Height 3" taller than selected Back Post Height.

## Back Upholstery

- 1920000 **No Back Upholstery**
- Standard Upholstery** Not available with Cane Handle Back. *E2611 / E2612* \$ 195
- 1920418 Short (18" Height)  1920421 Regular (21" Height)  1020424 Tall (24" Height)
- 1920450 **Tension Adjustable Upholstery** Height of Upholstery is 16". Not available with Cane Handle Back. *E2611 / E2612* \$ 295
- Position of Upholstery** Position of top of Upholstery relative to the Seat Pan. Circle below of select from the drop-down menu on the right.
- | Upholstery  | 22" Back Post |     |     |     |     |     |  |  |
|-------------|---------------|-----|-----|-----|-----|-----|--|--|
| Short       | 17"           | 18" | 19" | 20" | 21" | 22" |  |  |
| Regular     |               |     |     | 20" | 21" | 22" |  |  |
| Tension Adj |               |     | 19" | 20" | 21" | 22" |  |  |
- | Upholstery  | 25" Back Post |     |     |     |     |     |     |     |     |
|-------------|---------------|-----|-----|-----|-----|-----|-----|-----|-----|
| Short       | 17"           | 18" | 19" | 20" | 21" | 22" | 23" | 24" | 25" |
| Regular     |               |     |     | 20" | 21" | 22" | 23" | 24" | 25" |
| Tall        |               |     |     |     |     |     | 23" | 24" | 25" |
| Tension Adj |               |     | 19" | 20" | 21" | 22" | 23" | 24" | 25" |
- Height \_\_\_\_\_

## Tilt/Recline Options

- Attendant Tilt only** Select trigger mounting location below.
- 1923011 Right Back Post  1923012 Left Back Post
- Tilt Limit** Optional. If no tilt limit is selected, the full tilt range will be available.
- 1923220 Set at 20°  1923230 Set at 30°  1923240 Set at 40°
- 1923310 **Reclining Backrest (Manual)** \$ 895  
30° maximum recline. Not compatible with Backrest Extension, Cane Handle Back, Spring Back Suspension, Fold Down Back, Tension Adjustable Upholstery.
- 1923400 **Power Tilt** 250 lbs Weight Capacity. 45° tilt range. Includes 2 batteries for off-chair charging. *E1002* \$ 2,295  
Minimum Seat Width is 16". Power Tilt is compatible with either O2 Holder or IV Pole on the same chair, not both.
- Toggle Switch  1923421 Installed on Right side  1923422 Installed on Left side
- Button Switch  1923431 Installed on Right side  1923432 Installed on Left side
- Stealth Egg Switch  1923440 Not installed
- Optional for Egg Switch:  1923445 Swing Away Mounting Hardware Not compatible with all Headrest makes and models. *E1028* \$ 625
- Specify Headrest to be used:* \_\_\_\_\_

### Front Rigging

<input type="checkbox"/>	1925000	<b>No Front Rigging</b>	<input type="checkbox"/>	1925001	<b>Sector Blocks only</b>		
<b>Easy-Out™ Front Rigging</b> Effortless insertion and removal. No Hanger Bracket remains on frame. _____							
<input type="checkbox"/>	1925410	Fixed 70°					\$ 95
<input type="checkbox"/>	1925420	Universal Angle Adjustable (0°-110°)	Set at angle: _____ (use increments of 10°)				\$ 295
Adjust angle by loosening one bolt. Tool-free leg-length adjustment.							
<b>Swing Away Front Rigging</b> Legacy Design. Not compatible with Easy-Out.							
<input type="checkbox"/>	1925130	Swing Away 60°	<input type="checkbox"/>	1925140	Swing Away 70°		
<input type="checkbox"/>	1925200	Swing Away Elevating Legrests	Includes Calf Pads.			E0990	\$ 375
<input type="checkbox"/>	1925300	1" Offset Hangers for Swing Away					\$ 75
<input type="checkbox"/>	1925110	<b>Standard Plug In 70° Front Rigging</b>	Legacy Design. Not compatible with Easy-Out.				


### Footplates

<input type="checkbox"/>	1926010	<b>Composite Footplates</b>					
<input type="checkbox"/>	1926020	<b>Aluminum Angle Adjustable Footplates</b>				K0040	\$ 155
<input type="checkbox"/>	1926030	<b>One Piece Fixed Footplate</b>	Available on Plug In Front Rigging only. Requires selection of Composite or Aluminum Angle Adjustable.				\$ 350
<input type="checkbox"/>	1926040	<b>One Piece Flip-up Footplate</b>	Requires selection of Composite or Aluminum Angle Adjustable.				\$ 395
<input type="checkbox"/>	1926050	<b>Interlocking Angle Adjustable Footplates</b>	Available for Seat Widths up to 21" only.			K0040	\$ 195

### Footrest/Legrest Accessories

<b>Calf Pads</b> Not compatible with Gel Padding.							
<input type="checkbox"/>	1927012	For Easy-Out™ Front Rigging	Includes reinforced premium Calf Pad mount.				\$ 200
<input type="checkbox"/>	1927010	For Swing Away Front Rigging					\$ 135
<input type="checkbox"/>	1927020	<b>Gel Padding</b>	Not compatible with Calf Pad or Easy-Out Front Rigging.				\$ 125
<b>Gel Knee Pads</b>							
<input type="checkbox"/>	1927032	For Easy-Out™ Angle Adjustable Front Rigging					\$ 200
<input type="checkbox"/>	1927030	For Swing Away ELR Front Rigging					\$ 100
<input type="checkbox"/>	1927040	<b>Length Adjustable Heel Loops</b>				E0951	\$ 50
<input type="checkbox"/>	1927050	<b>Bodypoint Padded Calf Strap</b>				K0038	\$ 195
<b>Miller's Dynamic Footrest Mounts</b> Not compatible with Easy-Out Front Rigging.							
<input type="checkbox"/>	1927060	For Standard Plug In Front Rigging					\$ 445
<input type="checkbox"/>	1927065	For Swing Away 70° Front Rigging					\$ 575

### Front Rigging Extension Tubes

<b>Regular Extension Tubes</b>							
<input type="checkbox"/>	1928010	Short	14"-18" for Easy-Out Front Rigging	13"-17" for Swing Away and Standard Plug In Front Rigging			
<input type="checkbox"/>	1928020	Long	18"-24" for Easy-Out Front Rigging	15"-19" for Swing Away and Standard Plug In Front Rigging			\$ 75
<b>Extra Strength Extension Tubes (Solid)</b> Available on Plug In Front Rigging only. Not compatible with Easy-Out Front Rigging.							
<input type="checkbox"/>	1928110	Standard (13"-17")	<input type="checkbox"/>	1928120	Long (15"-19")		\$ 100
	<b>Lower Leg Length</b> If provided, Front Rigging will be set to this length by PDG. If left blank, the default setting will be used. Seat Pan to Foot Plates: _____ Measured along the length of the Front Rigging.						

### Miscellaneous Options

<input type="checkbox"/>	1930080	<b>Gel Hip Guides</b>	Reduces available Seat Width by 2".			E0956 & E1028	\$ 375
<input type="checkbox"/>	1930110	<b>IV Pole</b>	Not compatible with HD Anti Tipper, Attendant Wheel Locks.			E1029	\$ 295
<input type="checkbox"/>	1930120	<b>O2 Holder</b>	Not compatible with HD Anti Tipper, Attendant Wheel Locks.			E2208	\$ 395
<input type="checkbox"/>	1930170	<b>Adhesive Back Hook Strip, 2" x 30"</b>	For Cushion mounting.				\$ 10



# STELLAR GL

**Order Form - USA**  
Prices in U.S. Dollars. January 1, 2020

## On-Chair Seating & Positioning



Refer to the separate Seating & Positioning order form for PDG's on-chair options. Standard PDG discounts apply. \_\_\_\_\_

## Order Notes

HCPCS codes provided are based on code verification by PDAC or PDG's interpretation of existing code definitions.  
Use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.