

STELLAR

45° of tilt manual




- Ideal for hand propulsion with accessible rear wheels
- Features adjustable seat height, seat depth and back angle
- Knees remain low in tilt—wheel under desks, into vans and access other environments with ease
- Available in Heavy Duty (HD) extra-wide sizes
- Short wheelbase provides excellent maneuverability
- Transit-ready option
- ISO tested



Stellar HD

pega
MOBILITY



Nine standard frame colors to choose from and custom colors also available. Contact us for details.



Multi-adjustable



Adjustable push bar option



Easily adjustable axle position



Actuator adjusts for optimum size



HD option

features

| | | | |
|---------------------------------|------------|----------------------|-----------------------|
| Seat width (standard) | 14" – 22" | Overall width | 12" & 16" rear wheels |
| Seat width (HD) | 22" – 32" | | 20" & 24" rear wheels |
| Seat depth (adjustable) | 16" – 22" | | seat width + 9.5" |
| Seat height (adjustable) | 14" – 20" | | seat width + 10.5" |
| Back angle (adjustable) | 90° – 120° | | |
| Back height | 18" – 26" | | |
| Dynamic tilt | 45° | | |

Overall length
Weight capacity
Chair weight
Transit option available
Warranty

30" – 37" (without front rigging)
 250 lbs, HD to 450 lbs
 40 lbs (standard); 50 lbs (HD)
 ANSI/RESNA WC-19 (Standard frame only)
 PDG's exclusive TGIF warranty (Thank Goodness It Fits)
 Lifetime warranty on the frame.

PDG recommends all wheelchairs be fitted to the consumer under the guidance of a health care professional knowledgeable of the particular needs of the consumer and the limitations of the product. US Patent # 6 447064 B1