

## Hip Inflection Insert



## Padded Adductor Strap



Upper Weight Limit - 175 lbs!  
Maximum Height - 72"

Minimum Weight - 65 lbs.  
Minimum Height - 48"

Shoulder and Thigh Width Maxi-  
mums are Unlimited

Seat Depth - Two Sizes  
15" and 18"

Super low profile seat simplifies  
transfer of child into vehicle.

Designed and manufactured in the  
United States by  
Merritt Manufacturing, Inc.

### THE CHURCHILL

Is the first Belt Positioning Booster  
for children with Special Health Care  
Needs weighing 65-175lbs



## THE CHURCHILL™

A Belt Positioning Booster  
Designed for Children with  
Special Health Care Needs.



*Unique positioning  
features keep the child in  
position while using the  
vehicle lap and shoulder  
belt.*

# Positioning options Include...

## Vest

For complete torso support, this vest is fully adjustable, easy to put on, and features lower attachment straps to prevent ride up. Available in 2 sizes.



## Chest Wrap

Provides upper torso support.  
Available in 2 sizes.



## EZ-Install Head Rest



Prevents lateral movement.

## Stay-Put Pommel



## Hip Inflection Pommel

To Prevent slouching.



Use with the EZ-Up cap to restrict forward motion.