

# The Churchill



**A belt positioning backless booster  
designed for children and adults  
65-175 lbs  
with special health care needs**

The Churchill must be installed with LATCH (upper tether anchor and lower seat anchors) or with the Attachment Bracket kit. All Vehicles 2003 and newer are equipped with LATCH in at least two rear seating positions, so there should be no installation issues for these vehicles. Many 2000 - 2002 models also have LATCH (See your vehicle owner's manual to verify if LATCH is present).

**If your vehicle does not have latch there is a  
Attachment Bracket Kit available.**

# Special Features

- Super low profile base for easy transfer
- Shoulder, hip and thigh width unlimited
- Available seat base depth 15" or 18"
- Optional belt guides to properly position belt on thighs
- Adjustable Vest (Two Sizes Available)
- Hip/Leg straps to prevent slouching
- Foam Abductor
- Hip inflection/comfort pommel
- Latch attachment kit available
- Hip inflection wedge
- Headrest with lateral head support
- EZ-Up Cap™ (Three sizes available)



# Child transfer to Churchill Base

- Easy transfer of child to vehicle
- Super low profile base
- Swing away abductor
- Shoulder, hip and thigh width unlimited
- 15" or 18" seat depth available



# Churchill Belt Guides



- Optional belt guides to position vehicle seat belt low on upper thighs
- For larger child or adult belt guides can be moved to side if not needed – this will allow more hip room

# Churchill Vest /Shoulder Belt Guide



Easy to use shoulder belt guides to keep shoulder belt off the neck and positioned properly over the chest



Easy to install and adjust vest for upper body control  
(Two vest sizes available)



# Hip Inflection Wedge



Hip inflection wedge



# Churchill Pommels



- Standard pommel for keeping legs apart and in position



- Comfort Pommel to keep legs apart and position legs into a more comfortable position

# Leg Straps to Control Slouching



- Adjustable leg straps (standard)
- Leg straps keep child / adult hips in position – prevent slouching
- Can be used with or without the pommel



# A-ductor Strap

- Optional A-ductor strap to position and control leg movement
- A-ductor strap can be used with or without the pommel



# Churchill Head Rest



- Easy to install headrest mounts to vehicle head rest
- Provides lateral head support and sleeping support aid
- When vehicle does not have a headrest, the Churchill headrest attaches to seat mounting strap

# Churchill EZ-Up Caps



Available in three sizes

- EZ-Up Cap adjustable to fit the child's head
- EZ-Up Cap attaches to velcro located on headrest to keep head from falling forward

# Latch Bracket Attachment Kit



**The Attachment Bracket Kit is designed for use in vehicles that do not have LATCH.**

The Churchill must be installed with LATCH or with the Attachment Bracket kit. All Vehicles 2003 and newer are equipped with LATCH in at least two rear seating positions, so there should be no installation issues for these vehicles and the Attachment Bracket Kit will not be necessary. Many 2000 - 2002 models also have LATCH (See your vehicle owner's manual to verify if LATCH is present).

For vehicles without LATCH, the Attachment Bracket Kit can only be used on vehicle seats that you can slide a ruler, or similar object, from the back side of the vehicle seat through the crack in the vehicle seat (where the seat bottom and seat back come together) through to the front side of the seat. This can usually be done in mini-vans and SUVs. Generally, the Attachment Bracket Kit cannot be installed in the back seat of sedans. Installation of the Attachment Bracket Kit does require that a hole be drilled through the vehicle's floor board.