



2018 Churchill Retail Price List For Sapphire Blue Fabric Color Only

The Churchill-Pediatric Positioning
Vehicle Restraint System and Accessories

See Important Statement Below About Installing the Churchill****

Part Number	Description	Retail Price
2000CS15-S	The Churchill Standard* with Vest and 15 Inch Seat Depth (Sapphire Blue Cover) – a pediatric positioning vehicle restraint system (Weight range: 65 to 175 lbs.; Height Range: 48 to 72 inches). One Positioning Vest is included in price. Please choose vest size below.	\$1,095.00
2000CS18-S	The Churchill Standard* with Vest and 18 Inch Seat Depth (Sapphire Blue Cover) – a pediatric positioning vehicle restraint system (Weight range: 65 to 175 lbs.; Height Range: 48 to 72 inches). One Positioning Vest is included in price. Please choose vest size below.	\$1,095.00
2000CH15-S	The Churchill with Positioning Harness*** and 15 Inch Seat Depth (Sapphire Blue Cover) FOR USE WITH OUR ESCAPE PROOF PRODUCTS – a pediatric positioning vehicle restraint system (Weight range: 65 to 175 lbs.; Height Range: 48 to 72 inches). One Positioning Harness is included in price. NOTE: Add the escape proof Buckle Guard and Chest Clip Guard listed below for the proficient escape artist.	\$995.00
2000CH18-S	The Churchill with Positioning Harness*** and 18 Inch Seat Depth (Sapphire Blue Cover) FOR USE WITH OUR ESCAPE PROOF PRODUCTS – a pediatric positioning vehicle restraint system (Weight range: 65 to 175 lbs.; Height Range: 48 to 72 inches). One Positioning Harness is included in price. NOTE: Add the escape proof Buckle Guard and Chest Clip Guard listed below for the proficient escape artist.	\$995.00
2000CVS-S	Small Positioning Vest (Sapphire Blue) - for chest circumference between 25 and 32 inches and belly circumference less than 32 inches. **Price charged only when ordering vest separately.	\$275.00**
2000CVM-S	Medium Positioning Vest (Sapphire Blue) - for chest circumference between 30 and 37 inches and belly circumference less than 37 inches. **Price charged only when ordering vest separately.	\$275.00**
2000CVL-S	Large Positioning Vest (Sapphire Blue) - for chest circumference between 35 and 44 inches and belly circumference less than 44 inches. **Price charged only when ordering vest separately.	\$275.00**
2000BG	Buckle Guard – ONLY FOR USE WITH 2000CH15-S or 2000CH18-S	\$150.00
2000CCG-S	Chest Clip Guard – ONLY FOR USE WITH 2000CH15-S or 2000CH18-S	\$150.00
	Optional Accessories (next page)	

2000EUHR-S	EZ-Up Head Rest (Sapphire Blue) – Supports head to limit side to side motion. To limit forward motion, choose a cap size from below: The cap is designed to Velcro to the head rest.	\$155.00
See note to the right regarding the two cap options.	Note: There are two options for cap material. Sapphire Blue or Khaki. The Khaki option is the coolest option for hot temperatures because it is single layer.	
1000EUCS-K	EZ-Up Cap – Size Small (Khaki) - (Head Circumference: 17.5 – 19.5 inches)	\$89.00
1000EUCS-S	EZ-Up Cap – Size Small (Sapphire Blue) - (Head Circumference: 17.5 – 19.5 inches)	\$89.00
1000EUCM-K	EZ-Up Cap – Size Medium (Khaki) - (Head Circumference: 19.5 – 21.5 inches)	\$89.00
1000EUCM-S	EZ-Up Cap – Size Medium (Sapphire Blue) - (Head Circumference: 19.5 – 21.5 inches)	\$89.00
1000EUCL-K	EZ-Up Cap – Size Large (Khaki) - (Head Circumference: 21 – 23.25 inches)	\$89.00
1000EUCL-S	EZ-Up Cap – Size Large (Sapphire Blue) - (Head Circumference: 21 – 23.25 inches)	\$89.00
2000FA-S	Foam Abductor (Sapphire Blue)	\$105.00
2000HICP-S	Hip Inflection Comfort Pommel (Sapphire Blue)	\$145.00
2000HIW-S	Hip Inflection Wedge (Sapphire Blue)	\$105.00
2000IC-15	Quick-Change Incontinent Cover for 15 Inch Seat Depth Base	\$125.00
2000IC-18	Quick-Change Incontinent Cover for 18 Inch Seat Depth Base	\$125.00
2000ADD	Adductor Strap	\$75.00
2000AB	Attachment Bracket Kit (for Vehicles without LATCH) – See Important Note Below****	\$85.00

<p>* The “Churchill Standard” includes the following features:</p> <ul style="list-style-type: none"> • Base • Lower Anchor Attachment Straps • Torso Support System • Anti-Slouching Leg Strap Assembly • Foam Padding • Fabric Cover • Positioning Vest (Small, Medium, or Large) • Tether Extension Assembly • Belt Positioning Guides <p>** One Vest is included in price of “Standard Churchill with Vest.” Retail Price of additional vests is \$275.00 each.</p>	<p>***The “Churchill with Harness” includes the following features:</p> <ul style="list-style-type: none"> • Base • Lower Anchor Attachment Straps • Positioning Harness • Foam Padding • Fabric Cover • Tether Extension Assembly • Belt Positioning Guides
--	--

****** The Attachment Bracket Kit is designed for use in vehicles that do not have LATCH. The Churchill must be installed with LATCH or with the Attachment Bracket kit. All Vehicles 2003 and newer are equipped with LATCH in at least two rear seating positions, so there should be no installation issues for these vehicles and the Attachment Bracket will not be necessary. Many 2000 - 2002 models also have LATCH (See your vehicle owner’s manual to verify if LATCH is present).**

For vehicles without LATCH, the attachment bracket can only be used on vehicle seats that you can slide a ruler, or similar object, from the back side of the vehicle seat through the crack in the vehicle seat (where the seat bottom and

2018 Churchill Retail Price List for Sapphire Blue Fabric Color Only: Effective January, 2018

seat back come together) through to the front side of the seat. This can usually be done in mini-vans and SUVs. Generally, the Attachment Bracket cannot be installed in the back seat of sedans. Installation of the Attachment Bracket Kit does require that a hole be drilled through the vehicle's floor board.

Shipping charges – The cost of shipping "Best Way" will be added to Purchase Orders.



2018 Churchill Retail Price List For Total Camo Fabric Color Only

The Churchill-Pediatric Positioning
Vehicle Restraint System and Accessories

See Important Statement Below About Installing the Churchill****

Part Number	Description	Retail Price
2000CS15-TC	The Churchill Standard* with Vest and 15 Inch Seat Depth (Total Camo Cover) – a pediatric positioning vehicle restraint system (Weight range: 65 to 175 lbs.; Height Range: 48 to 72 inches). One Positioning Vest is included in price. Please choose vest size below.	\$1,095.00
2000CS18-TC	The Churchill Standard* with Vest and 18 Inch Seat Depth (Total Camo Cover) – a pediatric positioning vehicle restraint system (Weight range: 65 to 175 lbs.; Height Range: 48 to 72 inches). One Positioning Vest is included in price. Please choose vest size below.	\$1,095.00
2000CH15-TC	The Churchill with Positioning Harness*** and 15 Inch Seat Depth (Total Camo Cover) <u>FOR USE WITH OUR ESCAPE PROOF PRODUCTS</u> – a pediatric positioning vehicle restraint system (Weight range: 65 to 175 lbs.; Height Range: 48 to 72 inches). One Positioning Harness is included in price. NOTE: Add the escape proof Buckle Guard and Chest Clip Guard listed below for the proficient escape artist.	\$995.00
2000CH18-TC	The Churchill with Positioning Harness*** and 18 Inch Seat Depth (Total Camo Cover) <u>FOR USE WITH OUR ESCAPE PROOF PRODUCTS</u> – a pediatric positioning vehicle restraint system (Weight range: 65 to 175 lbs.; Height Range: 48 to 72 inches). One Positioning Harness is included in price. NOTE: Add the escape proof Buckle Guard and Chest Clip Guard listed below for the proficient escape artist.	\$995.00
2000CVS-TC	Small Positioning Vest (Total Camo) - for chest circumference between 25 and 32 inches and belly circumference less than 32 inches. **Price charged only when ordering vest separately.	\$275.00**
2000CVM-TC	Medium Positioning Vest (Total Camo) - for chest circumference between 30 and 37 inches and belly circumference less than 37 inches. **Price charged only when ordering vest separately.	\$275.00**
2000CVL-TC	Large Positioning Vest (Total Camo) - for chest circumference between 35 and 44 inches and belly circumference less than 44 inches. **Price charged only when ordering vest separately.	\$275.00**
2000BG	Buckle Guard – ONLY FOR USE WITH 2000CH15-TC or 2000CH18-TC	\$150.00
2000CCG-TC	Chest Clip Guard – ONLY FOR USE WITH 2000CH15-TC or 2000CH18-TC	\$150.00
	Optional Accessories (next page)	

2000EUHR-TC	EZ-Up Head Rest (Total Camo) – Supports head to limit side to side motion. To limit forward motion, choose a cap size from below: The cap is designed to Velcro to the head rest.	\$155.00
1000EUCS-TC	EZ-Up Cap – Size Small (Total Camo) - (Head Circumference: 17.5 – 19.5 inches)	\$89.00
1000EUCM-TC	EZ-Up Cap – Size Medium (Total Camo) - (Head Circumference: 19.5 – 21.5 inches)	\$89.00
1000EUCL-TC	EZ-Up Cap – Size Large (Total Camo) - (Head Circumference: 21 – 23.25 inches)	\$89.00
2000FA-TC	Foam Abductor (Total Camo)	\$105.00
2000HICP-TC	Hip Inflection Comfort Pommel (Total Camo)	\$145.00
2000HIW-TC	Hip Inflection Wedge (Total Camo)	\$105.00
2000IC-15	Quick-Change Incontinent Cover for 15 Inch Seat Depth Base	\$125.00
2000IC-18	Quick-Change Incontinent Cover for 18 Inch Seat Depth Base	\$125.00
2000ADD	Adductor Strap	\$75.00
2000AB	Attachment Bracket Kit (for Vehicles without LATCH) – See Important Note Below****	\$85.00

<p>* The "Churchill Standard" includes the following features:</p> <ul style="list-style-type: none"> • Base • Lower Anchor Attachment Straps • Torso Support System • Anti-Slouching Leg Strap Assembly • Foam Padding • Fabric Cover • Positioning Vest (Small, Medium, or Large) • Tether Extension Assembly • Belt Positioning Guides <p>** One Vest is included in price of "Standard Churchill with Vest." Retail Price of additional vests is \$275.00 each.</p>	<p>***The "Churchill with Harness" includes the following features:</p> <ul style="list-style-type: none"> • Base • Lower Anchor Attachment Straps • Positioning Harness • Foam Padding • Fabric Cover • Tether Extension Assembly • Belt Positioning Guides
--	--

**** The Attachment Bracket Kit is designed for use in vehicles that do not have LATCH. The Churchill must be installed with LATCH or with the Attachment Bracket kit. All Vehicles 2003 and newer are equipped with LATCH in at least two rear seating positions, so there should be no installation issues for these vehicles and the Attachment Bracket will not be necessary. Many 2000 - 2002 models also have LATCH (See your vehicle owner's manual to verify if LATCH is present).

For vehicles without LATCH, the attachment bracket can only be used on vehicle seats that you can slide a ruler, or similar object, from the back side of the vehicle seat through the crack in the vehicle seat (where the seat bottom and seat back come together) through to the front side of the seat. This can usually be done in mini-vans and SUVs. Generally, the Attachment Bracket cannot be installed in the back seat of sedans. Installation of the Attachment Bracket Kit does require that a hole be drilled through the vehicle's floor board.

Shipping charges – The cost of shipping "Best Way" will be added to Purchase Orders.



P.O. Number:

Date:

2018 Churchill Order Form Only For Total Camo Fabric Color

Billing Address

Contact Name:
Organization:
Street Address:
City, State:
ZIP Code:
Phone: Fax:
E-Mail:

Shipping Address (If Different)

Contact Name:
Organization:
Street Address:
City, State:
ZIP Code:
Phone: Fax:
E-Mail:

Qty	Part Number	Description	Unit Price	Line Total
	2000CS15-TC	The Churchill Standard* with Vest and 15 Inch Seat Depth (Total Camo Cover) – a pediatric positioning vehicle restraint system (Weight range: 65 to 175 lbs.; Height Range: 48 to 72 inches). One Positioning Vest is included in price. Please choose vest size below.		
	2000CS18-TC	The Churchill Standard* with Vest and 18 Inch Seat Depth (Total Camo Cover) – a pediatric positioning vehicle restraint system (Weight range: 65 to 175 lbs.; Height Range: 48 to 72 inches). One Positioning Vest is included in price. Please choose vest size below.		
	2000CH15-TC	The Churchill with Positioning Harness*** and 15 Inch Seat Depth (Total Camo Cover) FOR USE WITH OUR ESCAPE PROOF PRODUCTS – a pediatric positioning vehicle restraint system (Weight range: 65 to 175 lbs.; Height Range: 48 to 72 inches). One Positioning Harness is included in price. NOTE: Add the escape proof Buckle Guard and Chest Clip Guard listed below for the proficient escape artist.		
	2000CH18-TC	The Churchill with Positioning Harness*** and 18 Inch Seat Depth (Total Camo Cover) FOR USE WITH OUR ESCAPE PROOF PRODUCTS – a pediatric positioning vehicle restraint system (Weight range: 65 to 175 lbs.;		

		Height Range: 48 to 72 inches). One Positioning Harness is included in price. NOTE: Add the escape proof Buckle Guard and Chest Clip Guard listed below for the proficient escape artist.		
	2000CVS-TC	Small Positioning Vest (Total Camo) - for chest circumference between 25 and 32 inches and belly circumference less than 32 inches. **Price charged only when ordering vest separately.		
	2000CVM-TC	Medium Positioning Vest (Total Camo) - for chest circumference between 30 and 37 inches and belly circumference less than 37 inches. **Price charged only when ordering vest separately.		
	2000CVL-TC	Large Positioning Vest (Total Camo) - for chest circumference between 35 and 44 inches and belly circumference less than 44 inches. **Price charged only when ordering vest separately.		
	2000BG	Buckle Guard – ONLY FOR USE WITH 2000CH15-S or 2000CH18-S		
	2000CCG-TC	Chest Clip Guard – ONLY FOR USE WITH 2000CH15-S or 2000CH18-S		
		Optional Accessories		
	2000EUHR-TC	EZ-Up Head Rest (Total Camo) – Supports head to limit side to side motion. To limit forward motion, choose a cap size from below: The cap is designed to Velcro to the head rest.		
	See note to the right regarding the two cap options.	Note: There are two options for cap material. Total Camo or Khaki. The Khaki option is the coolest option for hot temperatures because it is single layer.		
	1000EUCS-K	EZ-Up Cap – Size Small (Khaki) - (Head Circumference: 17.5 – 19.5 inches)		
	1000EUCS-TC	EZ-Up Cap – Size Small (Total Camo) - (Head Circumference: 17.5 – 19.5 inches)		
	1000EUCM-K	EZ-Up Cap – Size Medium (Khaki) - (Head Circumference: 19.5 – 21.5 inches)		
	1000EUCM-TC	EZ-Up Cap – Size Medium (Total Camo) - (Head Circumference: 19.5 – 21.5 inches)		
	1000EUCL-K	EZ-Up Cap – Size Large (Khaki) - (Head Circumference: 21 – 23.25 inches)		
	1000EUCL-TC	EZ-Up Cap – Size Large (Total Camo) - (Head Circumference: 21 – 23.25 inches)		
	2000FA-TC	Foam Abductor (Total Camo)		
	2000HICP-TC	Hip Inflection Comfort Pommel (Total Camo)		
	2000HIW-TC	Hip Inflection Wedge (Total Camo)		

	2000IC-15-TC	Quick-Change Incontinent Cover for 15 Inch Seat Depth Base (Total Camo)		
	2000IC-18-TC	Quick-Change Incontinent Cover for 18 Inch Seat Depth Base (Total Camo)		
	2000ADD-TC	Adductor Strap (Total Camo)		
	2000AB	Attachment Bracket Kit (for Vehicles without LATCH) – See Important Note Below****		
			Sub Total**	
Indiana Residents Only; Not Applicable to other States			Sales Tax (7%)	
			TOTAL	

<p>* The "Churchill Standard" includes the following features:</p> <ul style="list-style-type: none"> • Base • Lower Anchor Attachment Straps • Position Harness for Vest Option • Anti-Slouching Leg Strap Assembly • **Positioning Vest (Small, Medium, or Large) • Foam Padding • Fabric Cover • Tether Extension Assembly • Belt Positioning Guides <p>** One Vest is included in price of "Standard Churchill with Vest." Retail Price of additional vests is \$275.00 each.</p>	<p>***The "Churchill with Harness" includes the following features:</p> <ul style="list-style-type: none"> • Base • Lower Anchor Attachment Straps • Positioning Harness for 5-Point Harness Option • 5-Point Harness • Foam Padding • Fabric Cover • Tether Extension Assembly • Belt Positioning Guides
--	--

**** The Attachment Bracket Kit is designed for use in vehicles that do not have LATCH. The Churchill must be installed with LATCH or with the Attachment Bracket kit. All Vehicles 2003 and newer are equipped with LATCH in at least two rear seating positions, so there should be no installation issues for these vehicles and the Attachment Bracket will not be necessary. Many 2000 - 2002 models also have LATCH (See your vehicle owner's manual to verify if LATCH is present).

For vehicles without LATCH, the attachment bracket can only be used on vehicle seats that you can slide a ruler, or similar object, from the back side of the vehicle seat through the crack in the vehicle seat (where the seat bottom and seat back come together) through to the front side of the seat. This can usually be done in mini-vans and SUVs. Generally, the Attachment Bracket cannot be installed in the back seat of sedans. Installation of the Attachment Bracket Kit does require that a hole be drilled through the vehicle's floor board.



P.O. Number:

Date:

2018 Churchill Order Form Only For Sapphire Blue Fabric Color

Billing Address

Contact Name:	
Organization:	
Street Address:	
City, State:	
ZIP Code:	
Phone:	Fax:
E-Mail:	

Shipping Address (If Different)

Contact Name:	
Organization:	
Street Address:	
City, State:	
ZIP Code:	
Phone:	Fax:
E-Mail:	

Qty	Part Number	Description	Unit Price	Line Total
	2000CS15-S	The Churchill Standard* with Vest and 15 Inch Seat Depth (Sapphire Blue Cover) – a pediatric positioning vehicle restraint system (Weight range: 65 to 175 lbs.; Height Range: 48 to 72 inches). One Positioning Vest is included in price. Please choose vest size below.		
	2000CS18-S	The Churchill Standard* with Vest and 18 Inch Seat Depth (Sapphire Blue Cover) – a pediatric positioning vehicle restraint system (Weight range: 65 to 175 lbs.; Height Range: 48 to 72 inches). One Positioning Vest is included in price. Please choose vest size below.		
	2000CH15-S	The Churchill with Positioning Harness*** and 15 Inch Seat Depth (Sapphire Blue Cover) FOR USE WITH OUR ESCAPE PROOF PRODUCTS – a pediatric positioning vehicle restraint system (Weight range: 65 to 175 lbs.; Height Range: 48 to 72 inches). One Positioning Harness is included in price. NOTE: Add the escape proof Buckle Guard and Chest Clip Guard listed below for the proficient escape artist.		
	2000CH18-S	The Churchill with Positioning Harness*** and 18 Inch Seat Depth (Sapphire Blue Cover) FOR USE WITH OUR ESCAPE PROOF PRODUCTS – a pediatric positioning vehicle restraint system (Weight range: 65 to 175 lbs.;		

		Height Range: 48 to 72 inches). One Positioning Harness is included in price. NOTE: Add the escape proof Buckle Guard and Chest Clip Guard listed below for the proficient escape artist.		
	2000CVS-S	Small Positioning Vest (Sapphire Blue) - for chest circumference between 25 and 32 inches and belly circumference less than 32 inches. **Price charged only when ordering vest separately.		
	2000CVM-S	Medium Positioning Vest (Sapphire Blue) - for chest circumference between 30 and 37 inches and belly circumference less than 37 inches. **Price charged only when ordering vest separately.		
	2000CVL-S	Large Positioning Vest (Sapphire Blue) - for chest circumference between 35 and 44 inches and belly circumference less than 44 inches. **Price charged only when ordering vest separately.		
	2000BG	Buckle Guard – ONLY FOR USE WITH 2000CH15-S or 2000CH18-S		
	2000CCG-S	Chest Clip Guard – ONLY FOR USE WITH 2000CH15-S or 2000CH18-S		
		Optional Accessories		
	2000EUHR-S	EZ-Up Head Rest (Sapphire Blue) – Supports head to limit side to side motion. To limit forward motion, choose a cap size from below: The cap is designed to Velcro to the head rest.		
	See note to the right regarding the two cap options.	Note: There are two options for cap material. Sapphire Blue or Khaki. The Khaki option is the coolest option for hot temperatures because it is single layer.		
	1000EUCS-K	EZ-Up Cap – Size Small (Khaki) - (Head Circumference: 17.5 – 19.5 inches)		
	1000EUCS-S	EZ-Up Cap – Size Small (Sapphire Blue) - (Head Circumference: 17.5 – 19.5 inches)		
	1000EUCM-K	EZ-Up Cap – Size Medium (Khaki) - (Head Circumference: 19.5 – 21.5 inches)		
	1000EUCM-S	EZ-Up Cap – Size Medium (Sapphire Blue) - (Head Circumference: 19.5 – 21.5 inches)		
	1000EUCL-K	EZ-Up Cap – Size Large (Khaki) - (Head Circumference: 21 – 23.25 inches)		
	1000EUCL-S	EZ-Up Cap – Size Large (Sapphire Blue) - (Head Circumference: 21 – 23.25 inches)		
	2000FA-S	Foam Abductor (Sapphire Blue)		
	2000HICP-S	Hip Inflection Comfort Pommel (Sapphire Blue)		
	2000HIW-S	Hip Inflection Wedge (Sapphire Blue)		

	2000C-15-S	Quick-Change Incontinent Cover for 15 Inch Seat Depth Base (Sapphire Blue)		
	2000C-18-S	Quick-Change Incontinent Cover for 18 Inch Seat Depth Base (Sapphire Blue)		
	2000ADD-S	Adductor Strap (Sapphire Blue)		
	2000AB	Attachment Bracket Kit (for Vehicles without LATCH) – See Important Note Below****		
			Sub Total**	
Indiana Residents Only; Not Applicable to other States			Sales Tax (7%)	
				TOTAL

<p>* The "Churchill Standard" includes the following features:</p> <ul style="list-style-type: none"> • Base • Lower Anchor Attachment Straps • Positioning Harness for Vest Option • Anti-Slouching Leg Strap Assembly • **Positioning Vest (Small, Medium, or Large) • Foam Padding • Fabric Cover • Tether Extension Assembly • Belt Positioning Guides <p>** One Vest is included in price of "Standard Churchill with Vest." Retail Price of additional vests is \$275.00 each.</p>	<p>***The "Churchill with Harness" includes the following features:</p> <ul style="list-style-type: none"> • Base • Lower Anchor Attachment Straps • Positioning Harness for 5-Point Harness Option • 5-Point Harness • Foam Padding • Fabric Cover • Tether Extension Assembly • Belt Positioning Guides
---	--

**** The Attachment Bracket Kit is designed for use in vehicles that do not have LATCH. The Churchill must be installed with LATCH or with the Attachment Bracket kit. All Vehicles 2003 and newer are equipped with LATCH in at least two rear seating positions, so there should be no installation issues for these vehicles and the Attachment Bracket will not be necessary. Many 2000 - 2002 models also have LATCH (See your vehicle owner's manual to verify if LATCH is present).

For vehicles without LATCH, the attachment bracket can only be used on vehicle seats that you can slide a ruler, or similar object, from the back side of the vehicle seat through the crack in the vehicle seat (where the seat bottom and seat back come together) through to the front side of the seat. This can usually be done in mini-vans and SUVs. Generally, the Attachment Bracket cannot be installed in the back seat of sedans. Installation of the Attachment Bracket Kit does require that a hole be drilled through the vehicle's floor board.