

Pedicraft Canopy Enclosed Bed Specifications

	<u>Regular Beds</u>	<u>Wide Beds</u>
Overall length:	76"	76"
Overall width including side rails:	40"	48"
Overall height:	75"	75"
With Leg Extensions:	85"	85"
Side Rail Height/Length	22"/62"	22"/62"
Distance:		
Top of mattress to floor:	27"	27"
With Leg Extensions:	37"	37"
Top of mattress to top of enclosure:	49"	49"
Mattress width:	33.5"	41.5"
Mattress thickness:	3" & 2"	3" & 2"
Mattress length 3":	77"	77"
With Padded Head/Foot:	71"	71"
Mattress length 2":	76"	76"
With Padded Head/Foot:	73"	73"

Head Elevation Option operates by hand crank. Highest head elevation angle: 30 degrees.

Height Adjustable models operate by a hydraulic pump system.

Low position: 27" from floor. High position: 37" from floor.

Patient weight capacity: 300 pounds

Foam padding covers all metal frame components. Sure-Chek® anti-microbial, fire retardant vinyl covers foam padding and is fastened with YKK locking zippers. High strength mesh is 100% polyester.

The mattress is constructed of high density fire-retardant foam, covered in Sure-Chek® Vinyl and complies with the Consumer Products Safety Commission standard for flammability, 16 CFR Part 1633.

Distance from the top of 3" mattress to the top of the side gate is 13".