

Pedicraft®

Reflux Wedges
and Slings



Pedicraft's Reflux Wedge

Gastro-esophageal reflux in infants is often treated simply by placing the patient at a 30° angle after feeding. But it's difficult to keep a baby in that position.

Pedicraft Reflux Wedges solve the problem in a way that's convenient for the hospital nursing staff, comfortable for Baby, and easy for Mom and Dad to use at home. You don't need a special crib with an inclined mattress pan, you don't even need to keep the baby in the crib—the Reflux Wedge can go wherever you go.



Pedicraft®

Keeping Babies & Children Safe
and Comfortable for over 30 Years



Reflux Wedge

Designed specifically for infants with gastro-esophageal reflux, the Reflux Wedge will hold the patient comfortably at a 30° angle. High density, fire retardant foam is covered with non absorbent Sure-Chek® vinyl. Quilted cover zips off easily for washing. Two removable 100% cotton flannel dribble pads come with each wedge. Additional dribble pads and wedge covers are available. Specify blue or pink.

- Pedicraft Wedges come in six sizes. Use the following chart as a guide to order the correct size to accommodate your needs.

<u>Size</u>	<u>Pounds</u>	<u>Kilograms</u>
Bassinet Wedge	up to 5	up to 3
Extra Small	up to 3.5	up to 2
Small	3.5–5.5	2–3
Medium	5–10	2.5–5.5
Large	9–22	5–10
Extra Large	20–30	10–17

Weight ranges are approximate.

Specify Code when ordering

#7358 Bassinet Wedge	#7362 Large Wedge
#7359 Extra Small Wedge	#7462 Large Wedge with Snaps
#7559 Extra Small Wedge Replacement Cover	#7562 Large Wedge Replacement Cover
#7360 Small Wedge	#7662 Large Wedge Replacement Cover with Snaps
#7560 Small Wedge Replacement Cover	#7364 Extra Large Wedge
#7361 Medium Wedge	#7563 Extra Large Wedge Replacement Cover
#7461 Medium Wedge with Snaps	#7341 One Dozen Dribble Pads (specify size)
#7561 Medium Wedge Replacement Cover	
#7661 Medium Wedge Replacement Cover with Snaps	



Reflux Sling

If you have the ability to incline your mattress to a 30° angle, the Reflux Sling will keep baby secure and keep him from sliding down in the crib. Made of quilted fabric for patient comfort and durability. The sling comes in pink or blue and ties easily to the crib bars. Machine washable.

#7379	Extra Small	up to 5 pounds	up to 3 kgs
#7380	Small	3 to 7 pounds	1.5 to 4 kgs
#7381	Medium	5 to 10 pounds	2.5 to 5 kgs
#7382	Large	9 to 22 pounds	5 to 10 kgs



Snaps at the edge of the diaper area are optional on Medium and Large Reflux Wedges. Standard on other sizes.

Note to Parents

The Pedicraft Reflux Wedge is a medical device. It is designed to hold a baby in the prone position (on his or her stomach) for the treatment of gastro-esophageal reflux. Your doctor may advise variations in treatment based on individual circumstances.

Follow your doctor's advice about the amount of time your baby needs to spend in the Wedge. The hour after each feeding is probably the most important.



Bassinet Pouch Wedge. #7358

Perfect for newborns to snuggle up in. Lower angle allows wedge to fit into isolette.

Reflux Wedge Homecare Package

Includes: Wedge, 2 Wedge Covers, 4 Dribble Pads

#7491	Medium
#7492	Large
#7493	Extra Large

