

NAME: _____

PHONE: _____

SCHOOL/ORGANIZATION: _____

E-MAIL: _____

DEALER/REP: _____

SALES REP: _____ QUOTE: _____



Custom build your aluminum treatment table to meet the specific needs of your practice. The extraordinary load capacity of 650 lbs and the welded aluminum frame make the Aluma Elite Basic Treatment Table the perfect fit for therapists and patients alike. Each table comes standard with a welded aluminum frame, H-brace construction, comfortable 2" high density foam cushion, and softly rounded corners. Make it your own by adding length or width, laminate shelves or a cabinet that can be mounted on either end. Patient-friendly customizations include a nose hole and 70° manual lift back.

STANDARD FEATURES:

- Load Capacity: 650 lbs.
- Welded Aluminum Frame
- H-brace Construction
- 2" High density foam cushion
- Rounded corners
- Tough powder coat finish

Aluma Elite Basic

TREATMENT TABLES



A4327LB-3078 with Shelf, Cabinet, Nose Hole, (2) Drawers, and Slip Cover

FLAT TOP TABLES:

- | | | | |
|---|---|---|---|
| <input type="checkbox"/> A4327-2772
27" Flat Cushion
27"W x 72"L x 31"H
Qty _____ | <input type="checkbox"/> A4327-3072
30" Flat Cushion
30"W x 72"L x 31"H
Qty _____ | <input type="checkbox"/> A4327-2778
27" Flat Cushion
27"W x 78"L x 31"H
Qty _____ | <input type="checkbox"/> A4327-3078
30" Flat Cushion
30"W x 78"L x 31"H
Qty _____ |
|---|---|---|---|

LIFT BACK TABLES:

- | | | | |
|--|--|--|--|
| <input type="checkbox"/> A4327LB-2772
27" Lift Back Cushion
27"W x 72"L x 31"H
Qty _____ | <input type="checkbox"/> A4327LB-3072
30" Lift Back Cushion
30"W x 72"L x 31"H
Qty _____ | <input type="checkbox"/> A4327LB-2778
27" Lift Back Cushion
27"W x 78"L x 31"H
Qty _____ | <input type="checkbox"/> A4327LB-3078
30" Lift Back Cushion
30"W x 78"L x 31"H
Qty _____ |
|--|--|--|--|

OPTIONS:

- AEB-BS** - Laminate Bottom Shelf
- AEB-SC** - Laminate Shelf and Open Cabinet
Cabinet Location:
 Head End Foot End
- AEB-SCD** - Laminate Shelf, Laminate Cabinet & (2) Laminate Drawers w/Black Wire Pull Handles
Cabinet Location:
 Head End Foot End
Drawer Location:
 Drawers on Right Drawers on Left
- AEB-NH** - Nose Hole in Cushion

Aluma Elite Basic

TREATMENT TABLES

STANDARD FRAME CHOICES:

- Light Gray Black

STANDARD SHELF & CABINET CHOICES:

- Fog Gray Black



A4327LB-3078 with Shelf, Cabinet, Nose Hole, and (2) Drawers

STANDARD UPHOLSTERY CHOICES:

- | | |
|--|--|
| <input type="checkbox"/> Royal Blue | <input type="checkbox"/> American Beauty |
| <input type="checkbox"/> Imperial Blue | <input type="checkbox"/> Lark |
| <input type="checkbox"/> Blue Ridge | <input type="checkbox"/> Burgundy |
| <input type="checkbox"/> Deep Violet | <input type="checkbox"/> Dove |
| <input type="checkbox"/> Forest | <input type="checkbox"/> Charcoal |
| <input type="checkbox"/> Black | <input type="checkbox"/> Mocha |
| <input type="checkbox"/> Mandarin Orange | |

CUSTOM BRANDING:

- Decal Custom Small**
Small Logo Decal
(6"W x 6"H) 
- Decal Custom Large**
Large Logo Decal
(12"W x 12"H) 

CUSTOM CUSHION BRANDING:

- A4451-27**
Logo Slip Cover
(Lift Back Only)
(27"W x 12"H) 
- A4451**
Logo Slip Cover
(Lift Back Only)
(30"W x 12"H)
- LOGO DB**
Deboss Logo
(12"W x 12"H)
DEBPLT12
Deboss Plate 

- LOGO CP**
Color Edge Print
(12"W x 12"H) 

COLOR EDGE PRINT UPHOLSTERY CHOICES:

- | | |
|---|------------------------------------|
| <input type="checkbox"/> Blue Ribbon | <input type="checkbox"/> Celestial |
| <input type="checkbox"/> Purple Passion | <input type="checkbox"/> Valentine |
| <input type="checkbox"/> Crimson | <input type="checkbox"/> Dolphin |
| <input type="checkbox"/> Black | |

****Color Edge Print upholstery colors do not match Standard Upholstery colors.**