

NAME: _____

PHONE: _____

SCHOOL/ORGANIZATION: _____

E-MAIL: _____

DEALER/REP: _____

SALES REP: _____ QUOTE: _____



ESSENTIAL THERA-P™ ELECTRIC TREATMENT TABLE

STANDARD FEATURES:

- Load Capacity: 500 lbs.
- 27"W x 77"L x 17-35"H
- Durable and stylish steel base (EPT9020/EPT9028)
- White powder coat base with single foot control pedal
- Whisper quiet lift system
- Height adjustment: 17-35"
- 27" Wide cushions in 7 different configuration options
- 2" High density foam cushions with rounded corners
- Individual locking casters
- Head section angle range: +40° to -90° (EPT77: +35° to -45°)
- Lower extremity angle range: 0 to 65° (EPT77: 0 to 75°)
- Armrest horizontal range: 0 to 45°
- Armrest vertical range: 0 to 7"
- 7 Section configuration includes electric elevating midsection with hand control
- ADA Compliant
- 2 Year warranty

NOTE: International (240V) option available



Ease of use, multiple options and a sleek but durable design make the Essential THERA-P line of electric treatment tables the affordable choice of therapists everywhere. Seven top options allow therapists to choose the table that best meets their needs and that ensures patients' most comfortable positioning. The Essential THERA-P electric treatment table is constructed of precision-cut flat and tube steel for durability and style. It also has a whisper-quiet lift system with elevating foot pedal. Designed with simplicity in mind, the Essential THERA-P electric treatment table offers a standard white powder coat base.

TABLE CONFIGURATIONS:

<input type="checkbox"/> EPT01 1 Section 27" Width	<input type="checkbox"/> EPT02 2 Section 27" Width	<input type="checkbox"/> EPT03 3 Section 27" Width	<input type="checkbox"/> EPT04 4 Section 27" Width	<input type="checkbox"/> EPT05 5 Section 27" Width	<input type="checkbox"/> EPT06 6 Section 27" Width	<input type="checkbox"/> EPT77 7 Section 27" Width
Qty _____	Qty _____	Qty _____	Qty _____	Qty _____	Qty _____	Qty _____

ESSENTIAL THERA-P™ ELECTRIC TREATMENT TABLE

TABLE FEATURES:



EPT01




EPT02



EPT03



EPT04



EPT05



EPT06



EPT77

OPTIONS:

- ☐ **EPT9021** - 500 lb Lifting Capacity (240V for EPT01, EPT02, EPT03, EPT04, EPT05, & EPT06)
- ☐ **EPT9029** - 500 lb Lifting Capacity (240V for EPT77)
- ☐ **PT303x** - 30"W Cushion Configuration
- ☐ **S1077** - Dual Foot Control Pedal
- ☐ **S1135** - Manual Elevating Midsection (Only available on EPT03 & EPT05)

- ☐ **EPTPH01** - Paper Roll Holder and Cutter (Only available on EPT01, EPT02, or EPT03)
- ☐ **S1078** - Dual Foot Bar Control
- ☐ **EPT10001** - Nose Hole Plug



Paper Roll Holder and Cutter (EPTPH01)

STANDARD UPHOLSTERY CHOICES:

- | | |
|--|--|
| <input type="checkbox"/> Royal Blue | <input type="checkbox"/> American Beauty |
| <input type="checkbox"/> Imperial Blue | <input type="checkbox"/> Lark |
| <input type="checkbox"/> Blue Ridge | <input type="checkbox"/> Burgundy |
| <input type="checkbox"/> Deep Phoenician | <input type="checkbox"/> Dove |
| <input type="checkbox"/> Forest | <input type="checkbox"/> Charcoal |
| <input type="checkbox"/> Black | <input type="checkbox"/> Mocha |

STANDARD BASE COLOR:

- ☐ White