

NAME: _____

PHONE: _____

SCHOOL/ORGANIZATION: _____

E-MAIL: _____

DEALER/REP: _____

SALES REP: _____ QUOTE: _____



Meet the specific needs of your practice with the many options available with the Essential Wood Treatment Table. Each table comes standard with a poplar hardwood frame, h-brace construction, comfortable 2" high density foam cushion, softly rounded corners and dual-bolted folding leg assembly for easy shipping. Make it your own by adding length, width, shelves or a cabinet that can be mounted on either end. Other options available include patient-friendly customizations consisting of a nose hole cushion and a 70° manual lift back with 6 position elevation.

STANDARD FEATURES:

- Load Capacity: 500 lbs.
- Solid poplar hardwood frame
- H-brace construction
- 2" High density foam cushion
- Rounded corners
- Dual-bolted folding leg assembly for convenient shipping

Essential Wood
TREATMENT TABLE



TEWT-3072 with Custom Stain

STANDARD TABLES:

FEWT-2772
27" Flat Cushion
27"W x 72"L x 31"H
Qty _____

FEWT-3072
30" Flat Cushion
30"W x 72"L x 31"H
Qty _____

FEWT-2778
27" Flat Cushion
27"W x 78"L x 31"H
Qty _____

FEWT-3078
30" Flat Cushion
30"W x 78"L x 31"H
Qty _____

LIFT BACK TABLES:

TEWT-2772
27" Lift Back Cushion
27"W x 72"L x 31"H
Qty _____

TEWT-3072
30" Lift Back Cushion
30"W x 72"L x 31"H
Qty _____

TEWT-2778
27" Lift Back Cushion
27"W x 78"L x 31"H
Qty _____

TEWT-3078
30" Lift Back Cushion
30"W x 78"L x 31"H
Qty _____

Essential Wood TREATMENT TABLE

OPTIONS:

- EWT Bottom Shelf** - Bottom Shelf
- EWT Cabinet** - Shelf and Open Cabinet
 - Head End** **Foot End**
- EWT Drawers** - Shelf, Cabinet & (2) Drawers with Black Wire Pull Handles
 - Head End** **Foot End**
- EWT-NH** - Nose Hole in Cushion
(Not Available with Paper Roll Holder)
- EPTPH01** - Paper Roll Holder
(On Foot End with Paper Cutter on Head End)



FEWT-2772 with Paper Roll Holder



TEWT-2772 with Paper Roll Holder

WOOD STAIN OPTION:

- EWT-Stain

WOOD STAIN CHOICES:

- Tungsten Black Beauty

STANDARD UPHOLSTERY CHOICES:

- | | |
|-------------------------------------|---------------------------------------|
| <input type="checkbox"/> Navy | <input type="checkbox"/> Hunter Green |
| <input type="checkbox"/> Agate | <input type="checkbox"/> Mushroom |
| <input type="checkbox"/> Burgundy | <input type="checkbox"/> Buff |
| <input type="checkbox"/> Purple | <input type="checkbox"/> Chocolate |
| <input type="checkbox"/> Teal | <input type="checkbox"/> Black |
| <input type="checkbox"/> Flint Grey | |

CUSTOM BRANDING:

- Decal Custom Small**
Small Logo Decal (6"W x 6"H) 
- Decal Custom Large**
Large Logo Decal (12"W x 12"H) 

CUSTOM CUSHION BRANDING:

- TIS 27 COVER/LOGO**
Logo Slip Cover (Lift Back Only)
(27"W x 12"H) 
- TIS 30 COVER/LOGO**
Logo Slip Cover (Lift Back Only)
(30"W x 12"H)

*Slip cover not available with paper roll holder



TEWT-2772 with Nose Hole, Bottom Shelf, Cabinet and Drawers on Foot End