

NAME: \_\_\_\_\_

PHONE: \_\_\_\_\_

SCHOOL/ORGANIZATION: \_\_\_\_\_

E-MAIL: \_\_\_\_\_

DEALER/REP: \_\_\_\_\_

SALES REP: \_\_\_\_\_ QUOTE: \_\_\_\_\_



# THERA-P™ ELECTRIC TREATMENT TABLE

## STANDARD FEATURES:

- Load Capacity: 650 lbs.
- 27"W x 77"L x 17-35"H
- Durable and stylish steel base (PT9020/PT9028)
- White powder coat base
- Whisper quiet lift system with single foot control pedal
- Height adjustment: 17-35"
- 27" Wide cushions in 7 different configuration options
- 2" High density foam cushions with rounded corners
- Low profile locking casters
- Head section angle range: +40° to -90° (PT1177: +35° to -45°)
- Lower extremity angle range: 0 to 65° (PT1177: 0 to 75°)
- Armrest horizontal range: 0 to 45°
- Armrest vertical range: 0 to 7"
- 7 Section configuration includes electric elevating midsection with hand control
- ADA Compliant

**NOTE:** International (240V) option available



PT1103

Ease of use, multiple options and a sleek but durable design make the THERA-P line of electric treatment tables the choice of therapists everywhere. Seven top options allow therapists to choose the table that best meets their needs and that ensures patients' most comfortable positioning. The THERA-P electric treatment table is constructed of precision-cut flat and tube steel for durability and style. It also has a whisper-quiet lift system (with elevating foot pedal) and easy-to-operate caster system for better ergonomics.

## TABLE CONFIGURATIONS:



**PT1101**  
1 Section  
27" Width

Qty \_\_\_\_\_



**PT1102**  
2 Section  
27" Width

Qty \_\_\_\_\_



**PT1103**  
3 Section  
27" Width

Qty \_\_\_\_\_



**PT1104**  
4 Section  
27" Width

Qty \_\_\_\_\_




**PT1105**  
5 Section  
27" Width

Qty \_\_\_\_\_



**PT1106**  
6 Section  
27" Width

Qty \_\_\_\_\_



**PT1177**  
7 Section  
27" Width

Qty \_\_\_\_\_

# THERA-P™ ELECTRIC TREATMENT TABLE

## TABLE FEATURES:



PT1101




PT1102



PT1103



PT1104



PT1105



PT1106



PT1177

## OPTIONS:


- PT9021** - 650 lb Lifting Capacity (240V for PT1101, PT1102, PT1103, PT1104, PT1105, & PT1106)
- PT9029** - 650 lb Lifting Capacity (240V for PT1177)
- PT303x** - 30"W Cushion Configuration
- S1077** - Dual Foot Control Pedal
- S1135** - Manual Elevating Midsection (Only available on PT1103 & PT1105)
- EPTPH01** - Paper Roll Holder and Cutter (Only available on PT1101, PT1102, or PT1103)
- S1078** - Dual Foot Bar Control
- PT10001** - Nose Hole Plug



Dual Foot Control Pedal (x2) (S1077)



Manual Elevating Midsection (S1135)




Dual Foot Bar Control (S1078)

















Paper Roll Holder and Cutter (EPTPH01)













## STANDARD BASE COLOR:

-  White

## CUSTOM BASE COLORS:

- |   |   |
|---|---|
| <input type="checkbox"/>  Black      | <input type="checkbox"/>  Green        |
| <input type="checkbox"/>  Light Gray | <input type="checkbox"/>  Forest Green |
| <input type="checkbox"/>  Charcoal   | <input type="checkbox"/>  Purple       |
| <input type="checkbox"/>  Royal Blue | <input type="checkbox"/>  Yellow       |
| <input type="checkbox"/>  Navy Blue  | <input type="checkbox"/>  Orange       |
| <input type="checkbox"/>  Red        | <input type="checkbox"/>  Brown        |
| <input type="checkbox"/>  Burgundy   | <input type="checkbox"/>  Bright Blue  |

## STANDARD UPHOLSTERY CHOICES:

- |  |  |
|--|--|
| <input type="checkbox"/>  Royal Blue      | <input type="checkbox"/>  American Beauty |
| <input type="checkbox"/>  Imperial Blue   | <input type="checkbox"/>  Lark            |
| <input type="checkbox"/>  Blue Ridge      | <input type="checkbox"/>  Burgundy        |
| <input type="checkbox"/>  Deep Phoenician | <input type="checkbox"/>  Dove            |
| <input type="checkbox"/>  Forest          | <input type="checkbox"/>  Charcoal        |
| <input type="checkbox"/>  Black           | <input type="checkbox"/>  Mocha           |

## BRANDING KITS:

### **BRANDING KIT 1 - S1170**

- Slip Cover
- Custom Table Base Color
- Custom Table Base Graphics

*\*Not available on PT1101, PT1104, or PT1106*

### **BRANDING KIT 2 - S1171**

- Slip Cover
- White Table Base Color
- Custom Table Base Graphics

*\*Not available on PT1101, PT1104, or PT1106*

### **BRANDING KIT 3 - S1173**

- Deboss Logo
- Custom Table Base Color
- Custom Table Base Graphics

*\*Not available on PT1104 or PT1106*

### **BRANDING KIT 4 - S1174**

- Deboss Logo
- White Table Base Color
- Custom Table Base Graphics

*\*Not available on PT1104 or PT1106*

### **BRANDING KIT 5 - S1176**

- Color Edge Print
- Custom Table Base Color
- Custom Table Base Graphics

*\*Not available on PT1104 or PT1106*

### **BRANDING KIT 6 - S1177**

- Color Edge Print
- White Table Base Color
- Custom Table Base Graphics

*\*Not available on PT1104 or PT1106*

## CUSTOM CUSHION BRANDING:

- PT3027**  
Logo Slip Cover  
(27"W x 12"H)



- PT3030**  
Logo Slip Cover  
(30"W x 12"H)








- LOGO DB**  
Deboss Logo  
(12"W x 12"H)  
**DEBPLT12**  
Deboss Plate



- LOGO CP**  
Color Edge Print  
(12"W x 12"H)



## COLOR EDGE PRINT UPHOLSTERY CHOICES:

- |   |  |
|---|--|
| <input type="checkbox"/>  Blue Ribbon    | <input type="checkbox"/>  Celestial |
| <input type="checkbox"/>  Purple Passion | <input type="checkbox"/>  Valentine |
| <input type="checkbox"/>  Crimson        | <input type="checkbox"/>  Dolphin   |
| <input type="checkbox"/>  Black          |  |

*\*\*Color Edge Print upholstery colors do not match Standard Upholstery colors.*

## CUSTOM BASE OPTIONS:

- S1178**  
Table Base Graphics
- PT9020-CC**  
Pick Custom Base Color
- S2000**  
Custom Table Base Powder Coat Color and Table Base Graphics



*PT1106 with Dual Footbars and Table Base Graphics*