

TelescoPik® & Mini®

Worlds First Adjustable Length Reacher

The latest innovation in design and patented technology by PikStik to bring you the most versatile reachers ever.

- Two Sizes: Adjustable Length Reaching Tools 30" to 44" (75cm to 110cm) or 20" to 26" (50cm to 66cm)
- Locking ring design makes length adjustment easy and secure.
- A strong aluminum shaft works indoors and out. Will Not Rust
- Rotating 5 1/2" (14 cm) Jaw Opening — Retrieve the biggest pieces from the smallest spaces.

MiniPik	T2026	20-26"
TelescoPik	T3044	30-44"



Locking Ring
Simple quarter turn of locking ring allows for easy adjustment.



Comfortable Grip
Ergonomic design makes pickup painless and easy!

