



POPULAS
furniture by AD·AS

KA-8546CR

Adjustable Therashape™ table with comfort curves



Height adjustable therapy table featuring multiple comfort curves with adjustability for individual positioning or multi-user environments.


FEATURES:

- Standard tabletop construction with durable laminate finish and vinyl edge
- Multiple comfort curves for enhanced accessibility
- Steel column with durable powder coat finish
- Custom designed extruded metal columns provide precision movement for fast, simple knob-adjustment
- 6 adjustable base levelers


TECHNICAL DATA:

- Standard table dimensions: 91" x 42"
- Comfort curve dimensions: 18"W x 10"D
- Vertical adjustment capability: 10"
- Standard base color: Almond, or Graphite Silver


MODEL	DIMENSIONS	ADJUSTMENT RANGE	START/END HEIGHT	LEG SPAN
KA-8546	91" X 42"	10"	23", 33"	30", 33", 30"



Top View



Front View



Side View