

EasyPulse POC

OXYGEN CONCENTRATOR



PrecisionMedical

The Difference is **Precision.**

EasyPulse POC

OXYGEN CONCENTRATOR

When you place a Precision Medical POC with a patient, you can rely on our durable, proven design to support your business.

Lightweight, durable and proven designs

3L and 5L pulse settings utilizing VPSA technology

With CMV, patients receive therapeutic oxygen in the first third of inspiration

Simple, one touch operation

Fully functional on AC or DC power

FAA approved for use on commercial airlines

Optional backpack making it comfortable for patients to wear all day long

Best in industry
5-year warranty
covering sieve bed
replacements

Lowest cost of
ownership

Battery packs with
up to 9.4 hours with
EasyPulse POC-5 and
11 hours for EasyPulse
POC-3



Precision Medical POCs

Precision Medical POCs can help your business by eliminating the high recurring costs of deliveries and the burden of managing cylinders. Other manufacturers use basic pressure only systems causing frequent sieve bed replacements. With our dependable POCs you can expect to minimize patient calls and enjoy lower cost of ownership.

Innovative VPSA Technology with CMV

Only Precision Medical POCs use Vacuum Pressure Swing Adsorption technology ensuring each sieve bed is perfectly regenerated every cycle. This results in clean sieve beds that produce high oxygen purity for a long period of time.

With controlled minute volume (CMV), your oxygen patients receive therapeutic oxygen in the first third of the inspiratory phase for maximum efficiency, regardless of breath rate.

Easy to Use

An easy to use push button operation allows patients to press one button to start their POC on the prescribed flow setting. Simple design reduces phone calls and troubleshooting.

Proven Battery Life

The internal and external batteries combine to provide up to 11 hours* of battery life. This duration allows patients to be mobile with the confidence that their oxygen needs will be met, at home as well as when traveling. CMV technology ensures battery duration stays consistent on all flow settings.

Recognized Design with Longest Product Life

EasyPulse POCs' product design is recognized by providers not only in the USA, but around the world for being the most dependable POC in their fleet. With the lowest return rate of any POC, trust the company that has been a leader in homecare and hospital respiratory care since 1982.

*Battery duration times are based on new, fully charged batteries used with the EasyPulse POC system. Battery duration times will decline with battery age and environmental operating conditions.



Ordering Information

EasyPulse Portable Oxygen Concentrator

3L: PM4130P • **5L:** PM4150P

includes: EasyPulse POC (3L or 5L) with AC and DC power supplies, cords, carry bag and cannula

EasyPulse POC with Backpack

3L: PM4130PK • **5L:** PM4150PK

includes: EasyPulse POC (3L or 5L) with AC and DC power supplies, cords, backpack, and cannula

EasyPulse POC with External Battery

3L: PM4130PEB • **5L:** PM4150PEB

includes: EasyPulse POC (3L or 5L) with AC and DC power supplies, cords, carry bag, cannula, and external battery

EasyPulse POC with Backpack and External Battery

3L: PM4130PEBK • **5L:** PM4150PEBK

includes: EasyPulse POC (3L or 5L) with AC and DC power supplies, cords, backpack, cannula, and external battery

Accessories

EasyPulse POC External Battery - 506760

EasyPulse POC AC Power Supply - 506629

EasyPulse POC AC Power Cord - 506961

EasyPulse POC DC Car Charger - 506750

EasyPulse POC Carry Bag

3L: 507947

5L: 506707

EasyPulse POC Backpack

3L: 508124

5L: 507303

EasyPulse POC Filter (Pack of 10) - 505847-10

EasyPulse POC Filter (Pack of 30) - 505847-30

Cannula - 504833

Product Specifications

Specifications are subject to change without prior notice.

	EasyPulse 3L (PM4130)	EasyPulse 5L (PM4150)
Oxygen concentrations	Up to 95.5%	Up to 95.5%
Sound level (on setting 2)	42 dBA	40 dBA
Height	8.5 in (21.6 cm)	10.1 in (25.7 cm)
Width	6.5 in (16.5 cm)	6.5 in (16.5 cm)
Depth	4.5 in (11.4 cm)	4.5 in (11.4 cm)
Weight (with internal battery)	4.9 lbs (2.22 kg)	6.6 lbs (2.98 kg)
Battery Duration (at any breath rate)	Setting 1: 5.5 hrs Setting 2: 3.8 hrs Setting 3: 2.8 hrs	Setting 1: 4.7 hrs Setting 2: 3.4 hrs Setting 3: 2.7 hours Setting 4: 2.0 hrs Setting 5: 1.6 hours
Battery Duration with optional external battery	Up to 11 hours	Up to 9.4 hours
Warranty	5-year warranty on concentrator and sieve beds 1-year warranty on battery	5-year warranty on concentrator and sieve beds 1-year warranty on battery