

Model	Go-Go® Sport
Part Number	3-wheel: S73 4-wheel: S74
Weight Capacity	325 lbs.
Maximum Speed ¹	Up to 4.7 mph
Ground Clearance ³	1.5" at motor / 0.75" at anti-tips
Turning Radius ³	3-wheel: 36" 4-wheel: 50.5"
Overall Length ³	3-wheel: 42" 4-wheel: 43"
Overall Width³	21.25"
Seat-to-Ground Height	20.5" - 22"
Front Tires	3-wheel: 3" x 8" solid 4-wheel: 3" x 9" solid
Rear Tires	3" x 9" solid
Range Per Charge ^{1,2}	3-wheel: Up to 15.9 miles 4-wheel: Up to 14.5 miles
Total Weight Without Batteries	3-wheel: 93 lbs. 4-wheel: 99 lbs.
Heaviest Piece When Disassembled ³	3-wheel: 35 lbs. (front section) 4-wheel: 44 lbs. (front section)
Standard Seating	Type: Compact, foldable Weight: 23 lbs. Material: Black vinyl Dimensions: Width: 17", Depth: 17"
Motor Output	24-volt DC motor, sealed transaxle
Brakes	Regenerative and electromechanical
Battery Requirements ⁴	(Two) 18 AH, 12 V, deep-cycle Weight: 14 lbs. each
Battery Charger	Off-board, 3.5 amp

HCPCS Code: K0801 (Group 1, Heavy Duty)

¹ Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This specification can be subject to a variance of (+ or -) 10%.

² Tested in accordance with ANSI/RESNA, MC Vol2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.

³ Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.

⁴ Battery weight may vary based on battery manufacturer.

PRIDE[®]
MOBILITY
LIVE YOUR BEST[®]

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

GOGOSPORT_061819

GOGO[®]
SPORT

**FEATURES...
PERFORMANCE...
...AND
TRAVEL READY!**



LIVE YOUR BEST[®]

pridemobility.com |    

FEATURES...PERFORMANCE... ...AND TRAVEL READY!



Compact, foldable seat with armrests

Feather-touch disassembly

Delta tiller with charger port

High-intensity LED headlight

Interchangeable shroud panels in red and blue



Two sets of easily changeable, red and blue colored panels.



Feather-touch disassembly for quick transport or storage.



Charger port conveniently located on side of tiller.

The Go-Go® Sport delivers high-performance and convenience on the go. With a 325 lb. weight capacity, a charger port in the tiller, standard front LED lighting and a maximum speed up to 4.7 mph., the Go-Go Sport is feature rich and travel ready.

Additional Features

- Convenient dual voltage charger can charge battery pack on-board or off-board
- Deck basket for additional storage (4-wheel only)
- Front frame-mounted seat post offers maximum stability
- Front tiller-mounted basket
- Long-lasting, high-intensity, front LED headlight offers optimal pathway illumination
- Microprocessor-based controller offers optimal power management and added safety features
- Pride's exclusive black, non-scutting tires

Colors (both color panel sets included)



Red



Blue

Accessories*

- Additional 18AH battery box
- Cane/crutch holder
- Cell phone holder
- Cup holder
- 50" Lap belt
- Oxygen tank holder
- Rear basket
- Rearview mirror
- Saddlebag
- Safety flag
- Under seat storage (assembly required)
- Walker holder
- Weather cover
- Wishbone crutch holder
- XLR USB charger
- Additional 18AH battery box

Military Patches



* Only one rear accessory can be used at a time.