

ACHIEVABLE POSITIONS



Lifted



Seated



Reading



Napping



Relaxing



True Infinite



Trendelenburg



Zero Gravity




Lay Flat



2-POSITION TV RECLINE

Lifted • Seated • Reading


The 2-position TV recline is achievable on all Pride® Power Lift Recliners. This is a comfortable position for performing many relaxing activities including reading, watching TV, working on your tablet or for conversation with family and friends.



3-POSITION FULL RECLINE

Lifted • Seated • Reading
Napping

The 3-position full recline is available on most Pride® Power Lift Recliner models and puts you in a position ideal for napping or relaxing.



TRUE-INFINITE POSITION

Lifted • Seated • Reading
Napping • Relaxing

The true-infinite position recline is perfect for taking a nap. Infinite-position Pride® Power Lift Recliners enable the backrest and footrest to move independently from one another so that a variety of positions can be reached.



TRENDELENBURG POSITION

The Trendelenburg position, available on the LC-525i, LC-580i and 515iL, elevates one's feet higher than his or her head for increased comfort and relaxation.



ZERO GRAVITY POSITION

Lifted • Seated • Reading
Napping • Relaxing
Trendelenburg • Zero gravity

The zero gravity position, which can be reached on infinite-position Pride® Power Lift Recliner models, positions one's thighs and lower legs higher than his or her torso.



INFINITE LAY FLAT POSITION

Can achieve infinite (two, three, infinite) and lay flat positions