

The Next Generation **Superstand** **HLT**

**Taking Standing
to New Heights!**

◆ **High Loading
Transfer**

◆ **Tool Free
Adjustment**

◆ **Soft Touch
Knobs**

◆ **Interchangeable
Masts**

◆ **Easy Standing
Positions**

◆ **Modular
Folding Frame**



Shown in supine

MULTI-POSITION

- ◆ **PRONE**
- ◆ **SUPINE**
- ◆ **UPRIGHT**



A Division Of Axiom Industries
All specifications subject to change without notice

Superstand HLT EI (Early Intervention)

Client Height: Supine 20"-38" Prone 22"-40"

Base: 19" wide x 28" long. Loading Height 27"

Client weight limit 65 pounds

Superstand HLT Size 1

Client Height: Supine 28"-48" Prone 30"-50"

Base: 21" wide x 37 3/4" long. Loading Height 28"

Client weight limit 125 pounds

Superstand HLT Size 2

Client Height: Supine 40"-60" Prone 40"-62"

Base: 25" wide x 45" long. Loading Height 29"

Client weight limit 150 pounds

The Art of Standing



Along with maintaining all of the superior positioning capabilities of the original Superstand, the Superstand HLT offers a horizontal, high-loading transfer surface to facilitate convenient and safe transfers of the client into the standing system.



- ◆ **Easy Loading**
- ◆ **Multiple Positions**
- ◆ **Growth and Portability**

The modular frame not only folds for easier transport and storage, but will easily configure to fit a client through most of their growth years.



A Division Of Axiom Industries
All specifications subject to change without notice