

SUPERSTAND

Pediatric Multi-Position Standing System

SUPERSTAND YOUTH

Client: _____ Age: _____
 Height: _____ Weight: _____ Dx: _____
 Address: _____
 City, St, Zip: _____
 Phone: _____
 Facility: _____
 Therapist: _____ Phone: _____

PO: _____ Date: _____
 Dealer: _____
 Contact: _____
 Address: _____
 City, St, Zip: _____
 Phone: _____ Fax: _____
 E-mail: _____

SUPERSTAND

Multi-Position Pediatric Standing System

SUPERSTAND YOUTH

Standard Features: Prone, Vertical and Supine capable; Hip and chest positioning pads with laterals; Headrest; Chinrest; Cotton support harness; Hands free standing angle adjustment; Width, height and depth adjustable knee system; Soft hip support pad; Individual planter-flexion, dorsi-flexion, rotation and height adjustable foot system; Adjustable heel loops; Toe straps; Adjustment lock out block and Base with 3" plastic casters.

<input type="checkbox"/> SUPERSTAND Client Height: Supine 26" - 48" Prone 30" - 50" \$2,995.00 Weight Limit 125 pounds Base: 21" wide x 40" deep		<input type="checkbox"/> SUPERSTAND YOUTH Client Height: 40" - 60" \$3,695.00 Weight Limit 150 pounds Base: 24" wide x 48" deep
---	---	---

Frame Color: Candy Red Candy Blue Candy Purple Neon Pink Neon Orange Shocker Green Razzberry White Black
 Upholstery Color: Red Blue Black Custom Upholstery: Customer Supplied Material Requested Color _____
 Harness Color: Black Denim Blue Denim Black Vinyl Blue Vinyl Red Vinyl

Lateral Pads and Brackets:

Superstand is standard with 3 1/2"H X 3 1/2"D Chest Laterals and 5 1/2"H X 4"D Hip Laterals with 90 Degree Brackets.
 Superstand Youth is standard with 6 1/2"H X 4 1/2"D Chest Laterals and 6 1/2"H X 4 1/2"D Hip Laterals with 90 Degree Brackets.

Chest Pad Brackets: 90° Std. 1" Offset 2" Offset Upgrade
 Hip Pad Brackets: 90° Std. 1" Offset 2" Offset Upgrade

Superstand	
90 Degree:	Width adjustment from 9" - 13"
1" Offset:	Width adjustment from 7" - 11"
2" Offset:	Width adjustment from 5" - 9"

Superstand Youth	
90 Degree:	Width adjustment from 11" - 15"
1" Offset:	Width adjustment from 9" - 13"
2" Offset:	Width adjustment from 7" - 11"

Options	<input type="checkbox"/> Upper Extremity Positioning Tray Assembly: Multi-Position tray with individual height, depth and angle adjustment.\$ 499.00
	<input type="checkbox"/> Tray Options: <input type="checkbox"/> Opaque tray overlay.....\$ 179.00 <input type="checkbox"/> Clear tray overlay.....\$ 179.00 <input type="checkbox"/> Tray easel.....\$ 249.00
	<input type="checkbox"/> Activity Bowl: Tray cut-out with 7" bowl.....\$ 195.00
	<input type="checkbox"/> Padded Tray Cover: Removable padded cover.....\$ 225.00
	<input type="checkbox"/> Size-Adjustable Knee Pads (pair): Allows knee pad width adjustment (from 2 1/2" - 4 1/2") for growth.....\$ 365.00
	<input type="checkbox"/> D-Ring Footstraps:\$ 75.00
	<input type="checkbox"/> Shoe Holders (pair): Formed ABS shoe holders with ankle and toe straps for foot positioning.....\$ 265.00 Size: <input type="checkbox"/> Sm (6" long x 4 1/4" wide) <input type="checkbox"/> Med (8" long x 4 7/8" wide) <input type="checkbox"/> Lrg (10 1/2" long x 6" wide)
	<input type="checkbox"/> Height Adjustable Pommel: Assists in maintaining lower extremity position.....\$ 225.00
	<input type="checkbox"/> Pommel W/ Rigid Hip Pad: Height adjustable pommel with height and depth adjustable rigid hip pad.....\$ 419.00
<input type="checkbox"/> Steel Caster Upgrade: High performance easy rolling steel casters with lock for both swivel and roll.....\$ 199.00	
<input type="checkbox"/> Padded Head Strap: Sheepskin-lined strap with Velcro closure.....\$ 99.00	

PRIME ENGINEERING

A Division of Axiom Industries Inc.
 All specifications subject to change without notice

SUPERSTAND

Pediatric Multi-Position Standing System

SUPERSTAND YOUTH



Custom Support Pads: Use this area only if you need custom-sized support pads on your Superstand. All custom pads must be quoted for price and manufacturing time.

Chest Pad: ___" Wide x ___" High Standard Pad is 10"w x 8"h for Superstand and 12"w x 15"h for Superstand Youth.


Hip Pad: ___" Wide x ___" High Standard Pad is 10"w x 7"h for Superstand and 12"w x 8"h for Superstand Youth.

Chest Laterals ___" Wide x ___" High Standard Pad is 3 1/2" d x 3 1/2" h for Superstand and 4 1/2" d x 6 1/2" h for Superstand Youth.


Hip Laterals ___" Wide x ___" High Standard Pad is 4" d x 5 1/2" h for Superstand and 4 1/2" d x 6 1/2" h for Superstand Youth.

Custom Tray Surface:

Tray: ___" Wide x ___" Deep Standard Tray is 22 1/2"w x 19"d for Superstand and 28"w x 24"d for Superstand Youth



Please supply any information for cutouts or modifications desired in this area or submit a paper or cardboard template. Additional information may be added in the special needs area below.



Client Measurements:

A: Thorax Width _____" B: Hip Width _____" C: Hip-to-Knee _____"
 D: Knee-to-Heel _____" E: Shoe Size _____" F: Shoulder Height _____"

Note: All Measurements Assume Symmetrical Standing and Sitting Posture.

Special Needs:

We are proud of our proven commitment to serving the needs of our customers with the finest standing products available since 1984. We encourage you to call us with your special needs, questions, or any item that we may assist with in your standing program. Our standing and positioning specialists will be happy to answer your questions or assist you with the design and fabrication of a unit that meets your needs. Please supply as much of the information in the client measurement area of this form as possible and list special needs in the following area to speed up the design process and / or response time.

PRIME ENGINEERING

A Division of Axiom Industries Inc.

All specifications subject to change without notice

Doc # 090118SSPLOF