

UP-RITE I, II, III

Pediatric Mid-Line Positioning Standing Frame

Client: _____ Age: _____
 Height: _____ Weight: _____ Dx: _____
 Address: _____
 City, St, Zip: _____
 Phone: _____
 Facility: _____
 Therapist: _____ Phone: _____

PO: _____ Date: _____
 Dealer: _____
 Contact: _____
 Address: _____
 City, St, Zip: _____
 Phone: _____ Fax: _____
 E-mail: _____

Standard Features: Height and depth adjustable clear table top; Adjustable heel block; Depth, height and width adjustable knee pads; One double and one single body support strap.


UP-RITE I
\$ 995.00



Client height up to 36"
 Base dimensions:
 20" wide x 20" long
 Tray dimensions:
 18" wide x 14 1/2" deep
 Knee pad cutout:
 3" wide x 2 3/4" deep
 Unit weight: 21lbs.

Up-Rite shown with optional caster base

UP-RITE II
\$ 1,295.00



Client height up to 48"
 Base dimensions:
 24" wide x 24" long
 Tray dimensions:
 21" wide x 17.5" deep
 Knee pad cutout:
 3" wide x 2 3/4" deep
 Unit weight: 27 lbs.

Up-Rite shown with optional caster base

UP-RITE III
\$ 1,395.00




Client height up to 60"
 Base dimensions:
 32" wide x 32" long
 Tray dimensions:
 24" wide x 21" deep
 Knee pad cutout:
 3 1/2" wide x 2 3/4" deep
 Unit weight: 44 lbs.

Up-Rite shown with optional caster base

- Options:**
- Wheeled caster base.....\$ 225.00
 - Size-adjustable knee pads (pair).....\$ 365.00
 - Double body support strap upgrade.....\$ 59.00
 - Chest laterals (pair).....\$ 199.00
 - Hip laterals (pair).....\$ 199.00
 - Knee straps (pair).....\$ 89.00

- Tray cutout with activity bowl.....\$ 195.00
- Tray overlay Clear Opaque.....\$ 179.00
- Foot straps (pair).....\$ 55.00
- D-Ring foot straps (pair).....\$ 75.00
- Adjustable head strap.....\$ 119.00
- Shoe holders (pair, small or large).....\$ 219.00



Client Measurements:

A: Thorax Width _____" B: Hip Width _____" C: Hip-to-Knee _____"
 D: Knee-to-Heel _____" E: Shoe Size _____" F: Shoulder Height _____"
 Note: All Measurements Assume Symmetrical Standing and Sitting Posture.