


Deluxe Hammock Sling



Multi-Purpose Seated Transfer Sling

- Provides support around the hips and sacral area and under the thighs up to the shoulders with optional head support
- Multiple loop and strap options can accommodate various seated and semi-reclined positions
- Divided leg supports that can be applied and removed in a seated or supine position
- Leg straps can be applied in various positions for different levels of support
- Built-in caregiver handles on each side add an additional level of control when moving and transferring user
- Small to medium commode opening

Materials



Mesh / Poly
Opt. Head Support ✓



PolySlip
Opt. Head Support ✓



Spacer
Opt. Head Support ✓

Applications




Up from the floor
*Requires Head Support



Higher up in bed
*Requires Head Support



From bed to wheelchair



From bed to commode



From wheelchair to toilet



Transfer to tub



Positioning in chair

Deluxe Hammock Sling

Part Numbers

Mesh / Poly

Part No.	Size	SWL
8H2700	Child	450 lb (204 kg)
8H2600	Junior	450 lb (204 kg)
8H2500	S	625 lb (283 kg)
8H2400	M	625 lb (283 kg)
8H2300	L	625 lb (283 kg)
8H2200	XL	1,000 lb (454 kg)

PolySlip

Part No.	Size	SWL
8H1700	Child	450 lb (204 kg)
8H1600	Junior	450 lb (204 kg)
8H1500	S	625 lb (283 kg)
8H1400	M	625 lb (283 kg)
8H1300	L	625 lb (283 kg)
8H1200	XL	1,000 lb (454 kg)

Spacer

Part No.	Size	SWL
8H4700	Child	625 lb (283 kg)
8H4600	Junior	625 lb (283 kg)
8H4500	S	625 lb (283 kg)
8H4400	M	625 lb (283 kg)
8H4300	L	625 lb (283 kg)
8H4200	XL	625 lb (283 kg)

Mesh / Poly with Head Support

Part No.	Size	SWL
8H2730	Child	450 lb (204 kg)
8H2630	Junior	450 lb (204 kg)
8H2530	S	625 lb (283 kg)
8H2430	M	625 lb (283 kg)
8H2330	L	625 lb (283 kg)
8H2230	XL	1,000 lb (454 kg)

PolySlip with Head Support

Part No.	Size	SWL
8H1720	Child	450 lb (204 kg)
8H1620	Junior	450 lb (204 kg)
8H1520	S	625 lb (283 kg)
8H1420	M	625 lb (283 kg)
8H1320	L	625 lb (283 kg)
8H1220	XL	1,000 lb (454 kg)

Spacer with Head Support

Part No.	Size	SWL
8H4730	Child	625 lb (283 kg)
8H4630	Junior	625 lb (283 kg)
8H4530	S	625 lb (283 kg)
8H4430	M	625 lb (283 kg)
8H4330	L	625 lb (283 kg)
8H4230	XL	625 lb (283 kg)

Mesh / Poly with Hip Strap

Part No.	Size	SWL
8H2701	Child	450 lb (204 kg)
8H2601	Junior	450 lb (204 kg)
8H2501	S	625 lb (283 kg)
8H2401	M	625 lb (283 kg)
8H2301	L	625 lb (283 kg)
8H2201	XL	1,000 lb (454 kg)

PolySlip with Hip Strap

Part No.	Size	SWL
8H1701	Child	450 lb (204 kg)
8H1601	Junior	450 lb (204 kg)
8H1501	S	625 lb (283 kg)
8H1401	M	625 lb (283 kg)
8H1301	L	625 lb (283 kg)
8H1201	XL	1,000 lb (454 kg)

Spacer with Hip Strap

Part No.	Size	SWL
8H4701	Child	625 lb (283 kg)
8H4601	Junior	625 lb (283 kg)
8H4501	S	625 lb (283 kg)
8H4401	M	625 lb (283 kg)
8H4301	L	625 lb (283 kg)
8H4201	XL	625 lb (283 kg)

Mesh / Poly with Head Support and Hip Strap

Part No.	Size	SWL
8H2731	Child	450 lb (204 kg)
8H2631	Junior	450 lb (204 kg)
8H2531	S	625 lb (283 kg)
8H2431	M	625 lb (283 kg)
8H2331	L	625 lb (283 kg)
8H2231	XL	1,000 lb (454 kg)

PolySlip with Head Support and Hip Strap

Part No.	Size	SWL
8H1721	Child	450 lb (204 kg)
8H1621	Junior	450 lb (204 kg)
8H1521	S	625 lb (283 kg)
8H1421	M	625 lb (283 kg)
8H1321	L	625 lb (283 kg)
8H1221	XL	1,000 lb (454 kg)

Spacer with Head Support and Hip Strap

Part No.	Size	SWL
8H4731	Child	625 lb (283 kg)
8H4631	Junior	625 lb (283 kg)
8H4531	S	625 lb (283 kg)
8H4431	M	625 lb (283 kg)
8H4331	L	625 lb (283 kg)
8H4231	XL	625 lb (283 kg)