


# Dual Access Sling



## Seated Transfer Sling for Toileting


- Provides support from under the thighs, behind the mid to low back, around the torso with a Velcro and belt enclosure
- Side retention straps are optional
- Divided, padded leg supports can be applied and removed in a seated position
- Large commode opening

## Materials



Poly

## Applications



From bed to commode



From wheelchair to toilet



Transfer to tub

## Part Numbers

### Poly

Part No.	Size	SWL
8B0700	Child	450 lb (204 kg)
8B0600	Junior	450 lb (204 kg)
8B0500	S	625 lb (283 kg)
8B0400	M	625 lb (283 kg)
8B0300	L	625 lb (283 kg)
8B0200	XL	1,000 lb (454 kg)

### Poly with Side Retention Straps

Part No.	Size	SWL
8C0700	Child	450 lb (204 kg)
8C0600	Junior	450 lb (204 kg)
8C0500	S	625 lb (283 kg)
8C0400	M	625 lb (283 kg)
8C0300	L	625 lb (283 kg)
8C0200	XL	1,000 lb (454 kg)