

Hammock Sling



Multi-Purpose Seated Transfer Sling

- Provides support around the hips and sacral area and under the thighs up to the shoulders with optional head support
- Multiple loop and strap options can accommodate various seated and semi-reclined positions
- Divided leg supports that can be applied and removed in a seated or supine position
- Leg straps can be applied in various positions for different levels of support
- Built-in caregiver handles on each side add an additional level of control when moving and transferring user
- Small commode opening

Materials



Mesh / Poly
Integrated Head Support✓



PolySlip
Integrated Head Support✓



Spacer
Integrated Head Support✓



Quilted
Integrated Head Support✓

Applications




Up from the floor
*Requires Head Support



Higher up in bed
*Requires Head Support



From bed to wheelchair



From bed to commode



From wheelchair to toilet



Transfer to tub



Positioning in chair