


# Universal Sling



## Multi-Purpose Seated Transfer Sling

- Provides support around the hips and sacral area and under the thighs up to the shoulders with optional head support
- Generous lateral trunk support
- Multiple loop and strap options can accommodate various seated and semi-reclined positions
- Divided leg supports can be applied and removed in a seated or supine position
- Leg straps can be applied in various positions for different levels of support
- Built-in caregiver handles on each side add an additional level of control when moving and transferring user
- Medium commode opening

## Materials



Mesh / Poly  
Opt. Head Support✓



Quilted  
Opt. Head Support✓



Mesh  
W. Head Support✓



PolySlip  
Opt. Head Support✓



Spacer  
Opt. Head Support✓



Padded  
Opt. Head Support✓



Disposable  
Integrated Head Support✓



Premium Disposable  
Integrated Head Support✓

## Applications




Up from the floor  
\*Requires Head Support



Higher up in bed  
\*Requires Head Support



From bed to wheelchair



From bed to commode



From wheelchair to toilet



Transfer to tub



Positioning in chair