


# Stand Aid Sling



## Sit to Stand Transfer Sling

- Designed to provide support for sit to stand transfers, toileting, and standing training
- Provides support from mid to low back, around the torso with velcro and buckle enclosure
- Optional non slip back available with velcro enclosure

## Materials



Poly

## Applications



Standing training



From bed to wheelchair



From bed to commode



From wheelchair to toilet



Rising – sitting down



Positioning in chair

## Part Numbers

### Poly

Part No.	Size	SWL
280070	Junior	600 lb (272 kg)
280071	S	600 lb (272 kg)
280072	M	600 lb (272 kg)
280073	L	600 lb (272 kg)
280074	XL	600 lb (272 kg)
280075	XXL	600 lb (272 kg)

### Non-Slip

Part No.	Size	SWL
521041	Junior	600 lb (272 kg)
521042	S	600 lb (272 kg)
521043	M	600 lb (272 kg)
521044	L	600 lb (272 kg)
521045	XL	600 lb (272 kg)
521046	XXL	600 lb (272 kg)