

# Protekt® Aire 9900

True Low Air Loss Mattress System with Alternating Pressure and Pulsation



HCPCS  
E0277

## FEATURES AND BENEFITS:

- 1275 liters per minute true low air loss.
- Provides effective true low air loss with alternating pressure and pulsation.
- State of the art micro-controller technology, monitors blower speed to set pressures.
- User-friendly control panel with dome keys, LED display for simple user interface. Displays patient comfort pressure level by weight.
- Alternating pressure cycles of 5, 10, 15, 20, 25 & 30 minutes allows clinicians to customize pressure redistribution and skin surface interface.
- Cell-on-cell mattress design prevents patient from bottoming out in case of power failure.
- Automatic fowler (seat inflate) recognizes when the patient is in the upright position and increases air flow in the middle bladders by 20% to prevent bottoming out.
- Audio/visual alarms in the event of power outage or low pressures. Internal memory always retains current settings during power outage, will resume desired operation upon power return.
- Lock out key completely locks the control panel, to eliminate tampering of the patient settings.
- Nylon top cover is fluid resistant and quilted for added comfort.
- Optional air bolster helps mitigate accidental roll outs.
- Optional raised side rails available on all sizes.
- Mattress dimensions: 36"x80"x8".
- Available in 42", 48", 54" & 60" widths.
- Weight capacity: 1,000 lbs.
- Warranty: 2 Year Non-Prorated.

81090-36	Protekt® Aire 9900
81090-36RR	Protekt® Aire 9900 with Raised Rails
81090-36AB	Protekt® Aire 9900 with Air Bolsters
81090-42	Protekt® Aire 9900 42"
81090-42AB	Protekt® Aire 9900 with Air Bolsters 42"
81090-48	Protekt® Aire 9900 48"
81090-48AB	Protekt® Aire 9900 with Air Bolsters 48"
81090-54	Protekt® Aire 9900 54"
81090-60	Protekt® Aire 9900 60"

 **Proactive**  
medical products