

# PROTEKT® DASH

## 32500 | STANDING TRANSFER

The Propekt® Dash is your best choice for fall prevention, reducing caregiver lifting injuries and promoting patient mobility. The Propekt® Dash offers ample weight capacity at 500 lbs. and includes an adjustable base which provides clearance around common obstacles such as recliners allowing access to patients.

- **500 lb. weight capacity.**
- **Dual shin pads for extra support.**
- **Ergonomically designed foot pedal offers hands free opening of base.**
- **Dual rear locking casters.**
- **Simple solution to compliment your electric stand assist lifters.**
- **First choice for any fall prevention program.**
- **Support strap available - attachment hooks for optional side support and back strap.**



### SPECIFICATION FOR ITEM # 32500

Weight Capacity	500 lbs.
Minimum Base Width	25.4"
Open Base Width	35"
Overall Height	43"
Minimum Seat Height	25"
Seat Width	20"
Overall Length	34"
Shin Pads	yes
Optional Support Strap	yes
Optional Seat Locks	yes

**WARRANTY:**  
Frame 3 years