

PROTEKT® 500 LIFT

33500 | POWER PATIENT LIFT

The high-performance Protekt® 500 Lift is the perfect blend of form and function for safe and easy patient transfers. This electric-powered, high quality patient transfer system provides an efficient method to safely move patients, preventing patient and caregiver injuries. Designed with the facility in mind, the lift is capable of floor lifts and fits underneath most beds in the lowest position, making it the perfect solution for transferring patients in small rooms with low beds.

- **500 lb. weight capacity.**
- **6 point spreader bar allows for multiple sling options, safer patient handling and a wider range of transfers.**
- **Capable of lifting from the floor.**
- **Ergonomically designed foot pedal offers hands free opening of base.**
- **Manual emergency lowering provides safe response to power loss.**
- **Dual rear locking casters.**



- **Clip free sling hookups improve safety and reduce downtime.**
- **Modular design allows easy parts replacement.**
- **Advanced electronics with auxiliary up/down lift controls.**

SPECIFICATION FOR ITEM # 33500

Weight Capacity	500 lbs.
Base Legs Length	45.4"
Base Closed External Width	25.6"
Base Open External Width	44"
Mast Height	48.5"
Base Height	4.7"
Min. Lifting Height	27.2"
Max. Lifting Height	70.9"
Net Weight	95 lbs.
Shipping Weight	110 lbs.

WARRANTY:
Electronics 2 years
Frame 3 years

PROTEKT® 600 LIFT

33600 | BARIATRIC POWER PATIENT LIFT

The high-performance Protekt® 600 Lift is the perfect blend of form and function for safe and easy patient transfers. This electric-powered, high quality patient transfer system provides an efficient method to safely move patients, preventing patient and caregiver injuries. Designed with the facility in mind, the lift is capable of floor lifts and fits underneath most beds in the lowest position, making it the perfect solution for transferring patients in small rooms with low beds.

- **600 lb. weight capacity.**
- **6 point spreader bar allows for multiple sling options, safer patient handling and a wider range of transfers.**
- **Capable of lifting from the floor.**
- **Ergonomically designed foot pedal offers hands free opening of base.**
- **Manual emergency lowering provides safe response to power loss.**
- **Dual rear locking casters.**

SPECIFICATION FOR ITEM # 33600

Maximum Capacity	600 lbs.
Base Legs Length	45.4"
Base Closed External Width	25.6"
Base Open External Width	44"
Mast Height	48.5"
Base Height	4.7"
Min. Lifting Height	27.2"
Max. Lifting Height	70.9"
Net Weight	99 lbs.
Shipping Weight	115 lbs.

WARRANTY:
Electronics 2 years
Frame 3 years



- **Clip free sling hookups improve safety and reduce downtime.**
- **Modular design allows easy parts replacement.**
- **Advanced electronics with auxiliary up/down lift controls.**