

# PROTEKT® 500 STAND

## 34500 | SIT-TO-STAND PATIENT LIFT

The Protekt® 500 lb. Sit-to-Stand patient lift has a higher weight capacity than most other competitors sit-to-stand lifts. With its 3" dual front casters it ensures safe and easy transfers. Our performance class actuator includes an emergency quick release.

HCPCS  
E0637

- 4 sling hooks allows for multiple sling options, safer patient handling and a wider range of transfers.
- Fast speed lift motor reduces time for the transfer.
- Battery lasts twice as long as competing brands for less downtime (24 volts).
- Ergonomically designed foot pedal offers hands free opening of base.
- 4" rear casters with brake & 3" dual front casters.
- Single piece molded kneepad with handlebar.
- Height & depth adjustable kneepad position.
- Full parts breakdowns and modular design on lift and electronics.
- Advanced electronics with auxiliary up/down lift controls.

**WARRANTY:**  
Electronics 2 years  
Frame 3 years  
Battery 1 year

### SPECIFICATION FOR ITEM # 34500

Weight Capacity	500 lbs.
Unit Weight	119 lbs.
Base Ext. Width - Closed	26.4"
Mast Height	41.7"
Base Length	41.3"
Front Caster Height	4.7"
Knee Pad Depth Range	2.8"
Knee Pad Height Range	20.2"-23.7"
Footplate Width	13.8"
Footplate Length	14.5"
Base Opened Turning Diameter	54.2"
Base Closed Turning Diameter	50.0"

Battery	LP-TC12-BPAK
Charger (Plug)	LP-TP7C-ADP
Charging Cradle (No plug)	LP-TBC2-CRDL
Hand Pendant	LP-TH10-PEN



# PROTEKT® 600 STAND

## 34600 | SIT-TO-STAND PATIENT LIFT

The Propekt® 600 lb. Sit-to-Stand has extra wide adjustable kneepads that will help meet a wider range of patient needs. To support the extra weight, the 600 Stand has a heavy duty boom to help lift and support the patient for smoother and safer transfers. The emergency quick release makes the 600 Stand the safest lift in any situation.

HCPCS  
E0637

**WARRANTY:**  
Electronics 2 years  
Frame 3 years  
Battery 1 year

### SPECIFICATION FOR ITEM # 34600

Weight Capacity	600 lbs.
Unit Weight	132 lbs.
Base Ext. Width - Closed	26.4"
Mast Height	41.7"
Base Length	41.3"
Front Caster Height	4.7"
Knee Pad Depth Range	2.8"
Knee Pad Height Range	20.2"-23.7"
Footplate Width	18.5"
Footplate Length	16.6"
Base Opened Turning Diameter	54.2"
Base Closed Turning Diameter	50.0"

Battery	LP-TC12-BPAK
Charger (Plug)	LP-TP7C-ADP
Charging Cradle (No plug)	LP-TBC2-CRDL
Hand Pendant	LP-TH10-PEN

- 4 sling hooks allows for multiple sling options, safer patient handling and a wider range of transfers.
- Fast speed lift motor reduces time for the transfer.
- Battery lasts twice as long as competing brands for less downtime (24 volts).
- Ergonomically designed foot pedal offers hands free opening of base.
- 4" rear casters with brake & 3" dual front casters.
- Single piece molded kneepad with handlebar.
- \* Height & depth adjustable kneepad position.
- Full parts breakdowns and modular design on lift and electronics.
- Advanced electronics with auxiliary up/down lift controls.

