

Prospect Designs

Alicia Prone Stander



The **Alicia Prone Stander** is designed to be easy to use and not take up a lot of space while providing comfortable and safe upright positioning for young children. The child is always at the lowest possible position while the stander is adjusted up or down to fit their height. The choice of adjustments and supports provide the child with the appropriate amount of help they need to experience weight bearing in a standing position.



FEATURES

Basic Stander includes:

- Two pair of lateral supports for trunk and hips
- Padded chest board
- Two 2" velcro straps
- Contoured knee supports
- Shoe holders
- Tray
- Locking casters


Available in three sizes: (APS1) (APS2) (APS3)

SPECIFICATIONS

	(APS1)	(APS2)	(APS3)
*-Footrest to top of chest board	22" to 29"	27" 37"	36" to 45"
-Width between laterals	6 1/2" to 9"	7 1/2" to 10 1/2"	8 1/2" to 12"
-Tray size	22" X 22"	22" X 23"	22" X 25"
-Base-width X length (fully extended)	24" X 32"	25" X 34"	30" X 45"

*Key Measurement

ACCESSORIES



Posterior Pelvic Support Gate
 Height and depth adjustable hinged gate provides firm posterior support and make transfers easier for caregivers. Three sizes correspond to the three sizes of standers (PPSG1) (PPSG2) (PPSG3)

Back Support Extension
 Flat pad that extends above the posterior pelvic gate and provides firm support to block the child from arching back (BKSE1)



Inset Activity Box
 Plastic tub inset into the tray can be used for a variety of activities. Includes clear plastic overlay. (ACTBX)

Back and Head Support Extension
 Flat pads for both shoulders and head. Provides even greater support to block trunk and head extension (BKSE2)

